NOTICE: Varian, Inc. was acquired by Agilent Technologies in May 2010. This document is provided as a courtesy but is no longer kept current and thus will contain historical references to Varian. For more information, go to www.agilent.com/chem.

Varian, Inc.
2700 Mitchell Drive
Walnut Creek, CA 94598-1675/usa

1200 LC/MS APCI Option

Atmospheric Pressure Chemical Ionization (APCI)

Operator's Manual

All rights reserved including the right of reproduction in whole or in part in any form. This document may be electronically reproduced, distributed, or printed in its entirety provided this copyright and statement are attached. Any modification or any other reproduction, distribution, or use of this document or portions hereof is strictly prohibited without the express written permission of Varian, Inc.

COPYRIGHT 2006 All rights reserved.

VARIAN

Quality Systems At Varian, Inc.

The ISO 9000 series standards were created in Geneva in 1987 to cut through a morass of conflicting quality definitions. These standards define a model for quality assurance systems in product design, development, manufacturing, installation, service, and customer support. They are now the worldwide quality assurance benchmark used to gauge the strength of a company's commitment to quality, and the value of its quality systems.

Various organizations around the world, such as the British Standards Institution (BSI), provide certified, objective auditors to scrutinize quality procedures, product development, manufacturing processes, and customer satisfaction programs. No company can claim ISO 9000 series registration unless it receives a stamp of approval from the demanding quality assessors of BSI or similar accredited examining body. ISO 9000 series registration constitutes an objective third-party report to determine the level of a supplier's commitment to quality.

In 1992, Varian, Inc., Analytical Instruments became registered to the most comprehensive of the ISO 9000 series standards — ISO 9001. ISO 9001 registration means that every stage of our quality system, including product development, manufacturing, final test, shipping, and parts and supplies has been rigorously examined against the most exacting set of internationally recognized standards. It means we live up to a standard of quality that you can count on today, and into the future. Our Quality System has received ISO 9001 certification number FM21797.

The quality systems that earned us ISO 9001 registration have direct benefits for our customers:

- ◆ We can speed instruments to you faster than ever before. Emergency orders can be processed even faster.
- ◆ We fill your orders promptly and completely.
- ◆ We have implemented a system of continuous feedback from our customers — we are aware of your needs today and tomorrow.
- ◆ We have improved your productivity by cutting systems failure rates in half and speeding service response time.
- ◆ We have embedded continuous improvement into the fabric of our organization so that we can achieve even higher levels of quality in the future.
- ◆ We are embedding GLP requirements into our products and services to help you meet your regulatory compliance requirements.

ISO 9001 registration is not enough. For us, quality is defined by our customers. We are not satisfied unless you are satisfied. We are striving to understand customer needs, using independent surveys, user groups, customer advisory boards, and our “Hallmark of Quality” response program, in addition to individual face-to-face customer contact. Our products and our processes are configured to meet those needs.

We know that you are seeking more than the most advanced processes and top-notch applications expertise. You want to join forces with a partner committed to delivering world-class quality, reliability, and value — on time, every time.

Our overriding aim is to be that partner.

Varian, Inc. Analytical Instrument Warranty

Hardware Products

All analytical instruments sold by Varian, Inc. are warranted to be free from defects in material and workmanship for the periods specified and in accordance with the terms on the face of Varian's quotation or as otherwise agreed upon in writing between Varian and the Customer. The warranty period begins on the date of **shipment** from Varian to the original Customer. However, where installation is paid for by the Customer or included in the purchase price, the warranty period begins upon completion of installation. If the Customer schedules **installation** to start later than 30 days after delivery or if such delay is caused through the Customer's inability to provide adequate facilities or utilities or through failure to comply with Varian's reasonable pre-installation instructions or through other omissions by Customer, then the warranty period starts on the 31st day from date of shipment. Moreover Varian will charge the Customer for labor and other expenses involved in making multiple or follow-up installation service calls.

Software Products

Where software is provided within the frame of a license agreement concluded between the Customer and Varian, any warranty shall be strictly in accordance with the terms of such agreement.

In the absence of a license agreement and unless an alternate warranty period is agreed upon in writing between Varian and the Customer, the warranty period is as specified on the face of Varian's quotation. Varian warrants such software products, if used with and properly installed on Varian hardware or other hardware as specified by Varian to perform as described in the accompanying Operator's Manual and to be substantially free of those defects which cause failure to execute respective programming instructions; however, Varian does not warrant uninterrupted or error-free operation.

Remedies

The sole and exclusive remedy under hardware warranty shall be **repair** of instrument malfunctions which in Varian's opinion are due or traceable to defects in original materials or workmanship or, at Varian's option, **replacement** of the respective defective parts, provided that Varian may as an alternative elect to **refund** an equitable portion of the purchase price of the instrument or accessory.

Repair or replacement under warranty does not extend the original warranty period.

Repair or replacement under warranty claims shall be made in Varian's sole discretion either by sending a Customer Support Representative to the site or by authorizing the Customer to return the defective accessory or instrument to Varian or to send it to a designated service facility. The Customer shall be responsible for loss or damage in transit and shall prepay shipping cost. Varian will return the accessory or instrument to the Customer prepaid and insured. Claims for loss or damage in transit shall be filed by the Customer. To correct software operation anomalies, Varian will issue software revisions where such revisions exist and where, in Varian's opinion, this is the most efficient remedy.

Limitation of Warranty

This **warranty does not cover** software supplied by the Customer, equipment and software warranted by another manufacturer or replacement of expendable items and those of limited life, such as but not limited to: Filters, glassware, instrument status lamps, source lamps, septa, columns, fuses, chart paper and ink, nebulizers, flow cells, pistons, seals, fittings, valves, burners, sample tubes, probe inserts, print heads, glass lined tubing, pipe and tube fittings, variable temperature dewars, transfer lines, flexible discs, magnetic tape cassettes, electron multipliers, filaments, vacuum gaskets, seats and all parts exposed to samples and mobile phases.

This **warranty shall be void** in the event of accident, abuse, alteration, misuse, neglect, breakage, improper operation or maintenance, unauthorized or improper modifications or tampering, use in an unsuitable physical environment, use with a marginal power supply or use with other inadequate facilities or utilities. Reasonable care must be used to avoid hazards.

This warranty is expressly in lieu of and excludes all other express or implied warranties, including but not limited to warranties of merchantability and of fitness for particular purpose, use or application, and all other obligations or liabilities on the part of Varian, unless such other warranties, obligations or liabilities are expressly agreed to in writing by Varian.

Limitation of Remedies and Liability

The remedies provided herein are the sole and exclusive remedies of the Customer. In no case will Varian be liable for incidental or consequential damages, loss of use, loss of production or any other loss incurred.

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Safety Information

Operating Instructions

This instruction manual is provided to help you establish operating conditions which will permit safe and efficient use of your equipment. Special considerations and precautions are also described in the manual, which appear in the form of **NOTES**, **CAUTIONS**, and **WARNINGS** as described below. It is important that you operate your equipment in accordance with this instruction manual and any additional information which may be provided by Varian. Address any questions regarding the safe and proper use of your equipment to your local Varian office.

NOTE

Information to aid you in obtaining optimal performance from your instrument.

CAUTION

Alerts you to situations that may cause moderate injury and/or equipment damage, and how to avoid these situations.

WARNING

Alerts you to potentially hazardous situations that could result in serious injury, and how to avoid these situations.

Warning Symbol

WARNING:
SHOCK HAZARD

WARNING:
CHEMICAL HAZARD

WARNING:
BURN HAZARD

WARNING:
EYE HAZARD

WARNING:
FIRE HAZARD

WARNING:
EXPLOSION HAZARD

WARNING:
RADIATION SOURCE

WARNING:
MOVING PARTS

Warning Description

Hazardous voltages are present inside instrument. Disconnect from main power before removing screw-attached panels.

Hazardous chemicals may be present. Avoid contact, especially when replenishing reservoirs. Use proper eye and skin protection.

Very hot or cryogenically cold surfaces may be exposed. Use proper skin protection.

Eye damage could occur either from flying particles, chemicals, or UV radiation. Use proper eye and face protection.

The potential for fire may be present. Follow manual instructions for safe operation.

The potential for explosion may exist because of type of gas or liquid used.

Ionizing radiation source is present. Follow manual instructions for safe operation.

Keep hands and fingers away.

VARIAN

General Safety Precautions

Follow these safety practices to ensure safe equipment operation.

- Perform periodic leak checks on all supply lines and pneumatic plumbing.
- Do not allow gas lines to become kinked or punctured. Place lines away from foot traffic and extreme heat or cold.
- Store organic solvents in fireproof, vented and clearly labeled cabinets so they are easily identified as toxic and/or flammable materials.
- Do not accumulate waste solvents. Dispose of such materials through a regulated disposal program and not through municipal sewage lines.

NOTICE: This instrument has been tested per applicable requirements of EMC Directive as required to carry the European Union CE Mark. As such, this equipment may be susceptible to radiation/interference levels or frequencies which are not within the tested limits.

WARNING

This instrument is designed for chromatographic analysis of appropriately prepared samples. It must be operated using appropriate gases and/or solvents and within specified maximum ranges for pressure, flows, and temperatures as described in this manual. If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

WARNING

It is the responsibility of the Customer to inform Varian Customer Support Representatives if the instrument has been used for the analysis of hazardous biological, radioactive, or toxic samples, prior to any instrument service being performed or when an instrument is being returned to the Service Center for repair.

Electrical Hazards

- Disconnect the instrument from all power sources before removing protective panels to avoid exposure to potentially dangerous voltages.
- When it is necessary to use a non-original power cord plug, make sure the replacement cord adheres to the color coding and polarity described in the manual and all local building safety codes.
- Replace blown fuses with fuses of the size and rating stipulated on the fuse panel or in the manual.
- Replace faulty or frayed power cords immediately with the same type and rating.
- Make sure that voltage sources and line voltage match the value for which the instrument is wired.

Compressed Gas Cylinders

- Store and handle compressed gases carefully and in strict adherence to safety codes.
- Secure cylinders to an immovable structure or wall.
- Store and move cylinders in an upright, vertical position. Before transport, remove regulators and install cylinder cap.
- Store cylinders in a well-ventilated area away from heat, direct sunshine, freezing temperatures, and ignition sources.
- Mark cylinders clearly so there is no doubt as to their contents.
- Use only approved regulators and connections.
- Use only connector tubing that is chromatographically clean (Varian Part Number 03-918326-00) and has a pressure rating significantly greater than the highest outlet pressure from the regulator.

GC Safety Practices

Exhaust System

No special exhaust ducting is necessary for GC detectors installed in a well-ventilated room except when the detectors are used to test hazardous chemicals. If you do install ducting:

- Use only fireproof ducting.
- Install a blower at the duct outlet.
- Locate duct intakes such that their vibration or air movement does not effect detector operation.
- Check periodically for proper operation of the duct.
- Ensure proper ventilation in lab area.

Radioactive Source Detectors

- Read carefully and comply with all NOTES, CAUTIONS, and WARNINGS in the Ni⁶³ ECD manual.
- Perform the tests for removable radioactive contamination described in the Ni⁶³ ECD manual.
- Comply with leak test schedules and procedures.

Burn Hazard

Heated or cryogenically cooled zones of gas chromatographs can remain hot or cold for a considerable time after instrument power is turned off. To prevent painful burns, ensure that all heated or cooled areas have returned to room temperature or wear adequate hand protection before you touch potentially hot or cold surfaces.

LC Safety Practices

High Pressure Hazard

- If a line ruptures, a relief device opens, or a valve opens accidentally under pressure, potentially hazardous high liquid pressures can be generated by the pump causing a high velocity stream of volatile and/or toxic liquids.
- Wear face protection when you inject samples or perform routine maintenance.
- Never open a solvent line or valve under pressure. Stop the pump first and let the pressure drop to zero.
- Use shatter-proof reservoirs capable of operating at 50-60 psi.
- Keep the reservoir enclosure closed when the reservoir is under pressure.
- Read and adhere to all NOTES, CAUTIONS, and WARNINGS in the manual.

Flash Chromatography

The operator should be familiar with the physico-chemical properties of the components of the mobile phase.

Keep solvents from direct contact with the polyurethane supply tubing as certain solvents will cause weakening and leaks with possible bursting.

All components of the system should be connected to a common power supply and common ground. This ground must be a true ground rather than a floating ground.

Non-polar solvents can develop a static charge when pumped through the system. All vessels that contain mobile phase (including tubing and collection vessels) must be grounded to dissipate static electricity.

Employ static measuring and static discharge devices (e.g., air ionizers) to safeguard against the buildup of static electricity.

Ultraviolet Radiation

Liquid chromatograph detectors that use an ultraviolet light source have shielding to prevent radiation exposure to personnel.

For continued protection:

- Ensure that protective lamp covers of variable and fixed wavelength detectors are in place during operation.
- Do not look directly into detector fluid cells or at the UV light source. When inspecting the light source or fluid cell, always use protective eye covering such as borosilicate glass or polystyrene.

The following is a Federal Communications Commission advisory: This equipment has been tested and found to comply with the limits of a Class A computing device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Spare Parts Availability

It is the policy of Varian to provide operational spare parts for any instrument and major accessory for a period of five (5) years after shipment of the final production run of that instrument. Spare parts will be available after this five (5) year period but on an *as available* basis. Operational spare parts are defined as those individual electrical or mechanical parts that are susceptible to failure during their normal operation. Examples include relays, lamps, temperature probes, detector elements, motors, etc. Sheet metal parts, structural members or assemblies and castings, printed circuit boards, and functional modules are normally capable of being rebuilt to like-new condition throughout their useful life and therefore will be supplied only on an *as available* basis after the final production run of the instrument.

Service Availability

Varian provides a variety of services to support its customers after warranty expiration. Repair service can be provided by attractively priced service contracts or on a time and material basis. Technical support and training can be provided by qualified personnel on both a contractual or as-needed basis.

Varian, Inc. Analytical Instruments Sales Offices

For Sales or Service assistance and to order Parts and Supplies, contact your local Varian office.

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www.varianinc.com

Contents

Overview of APCI.....	3
Introduction	3
Basic Principles	3
Positive Ion APCI	4
Negative Ion APCI.....	4
Applications	5
Installation of APCI Chamber	6
APCI Daily Startup and Shutdown.....	8
Startup.....	8
Shutdown	8
Tuning	8
APCI Mobile Phase Selection.....	8
Optimizing S/N for APCI	9
Typical APCI Operating Parameters	9
APCI Performance Evaluation.....	10
Operating Conditions for Testing the APCI.....	10
APCI Sensitivity Test.....	10
 APCI General Maintenance	 11
Cleaning.....	11
Cleaning the Spray Shield.....	11
Cleaning the Vacuum Capillary.....	11
Cleaning the Corona Needle.....	12
Replacing and Aligning the Corona Needle	12
Replacing the Nebulizing Needle	13
Removal	13
Installation	13
Replacing the APCI Quartz Tube	13
Removing the Quartz Tube	14
Installing the Quartz Tube	15
APCI Troubleshooting	17
Comparison of ESI and APCI Parameters	17
 APCI Spare Parts	 19

Overview of APCI

Introduction

For many small molecules, atmospheric pressure chemical ionization (APCI) is an excellent alternative to electrospray ionization (ESI). APCI assumes the analyte molecule has reasonable thermal stability. Thermal degradation in APCI is less likely than in GC/MS, but thermal stability is a fundamental restriction for APCI that is not shared by ESI. Optimization of APCI also requires an understanding of the relationship between the mobile phase solvents and the chemistry in the ionization process. Unlike ESI in which ions are extracted from droplets of mobile phase, APCI is a gas phase reaction between a reagent ion (typically from the mobile phase) and a neutral molecule of the analyte. This difference alters the APCI conditions for reversed phase chromatography and opens APCI to other modes of chromatography such as normal phase, and, in turn, to modified sample preparation schemes.

APCI and ESI are complementary, but significantly different in terms of setup and optimization. When the appropriate criteria for APCI have been met, chemical ionization offers advantages of increased sensitivity, enhanced selectivity, reduced interference, and robust, precise performance. For a significant number of small molecule applications, APCI may be the better solution than ESI.

Basic Principles

Atmospheric Pressure Chemical Ionization (APCI) requires the following:

- The liquid flow is rapidly vaporized
- solvent molecules are ionized in the corona discharge
- Chemical Ionization ionizes sample molecules

The APCI process begins with pneumatically assisted nebulization of the LC eluent. The spray of eluent droplets are completely vaporized in a hot, turbulent flow of auxiliary gas. The stream of vaporized mobile phase and sample components flows into the API chamber and crosses the path of a corona discharge. Through a multi-step process, the corona discharge creates CI reagent ions from the large excess of mobile phase solvent molecules. A final ion-molecule reaction between these CI reagent ions and neutral molecules of analyte generates the analyte ions the mass spectrometer monitors.

Positive Ion APCI

Positive Ion APCI can occur as follows:


- Proton Transfer: $\text{H}_3\text{O}^+ + \text{M} \rightarrow (\text{M}+\text{H})^+ + \text{H}_2\text{O}$
 - Also CH_3OH_2^+ , CH_3CNH^+ , NH_4^+
- Adduct Attachment: $\text{NH}_4^+ + \text{M} \rightarrow (\text{M} + \text{NH}_4)^+$
 - Also CH_3OH_2^+ , CH_3CNH^+ , H_3O^+
- Charge Exchange

Negative Ion APCI

Negative Ion APCI can occur as follows:

- Proton Abstraction: $\text{OH}^- + \text{M} \rightarrow (\text{M} - \text{H})^- + \text{H}_2\text{O}$
 - Also CH_3O^- , CH_2CN^-
- Adduct Attachment: $\text{CH}_3\text{COO}^- + \text{M} \rightarrow (\text{M} + \text{CH}_3\text{COO})^-$
 - Also Cl^- , HCOO^-
- Electron Capture

The APCI Source consists of two distinct parts, the atmospheric pressure chamber and the API capillary assembly. Vaporization, ionization and desolvation occur in the atmospheric chamber. The API capillary assembly uses a counter flow of drying gas to decluster heavily solvated ions. Using optimized gas flows and electrical potential, the API interface admits ions from the CI reaction while preventing neutrals from entering the sampling capillary.

Cutaway View of APCI Chamber

APCI is a soft ionization process that typically produces only singly charged ions with information about the molecular ion. Adduct ions are common, but lower mass fragment ions are not. The intense pseudo-molecular ion is, however, an ideal precursor ion for MS/MS. The requirements for volatility and thermal stability limit application of APCI to molecules below ~1500 daltons.

Applications

APCI is ideally suited to non-polar and medium polarity molecules. Analytes such as polynuclear aromatic hydrocarbons (PAHs), polychlorinated biphenyls (PCBs), fatty acids, phthalates, and triglycerides are perfect candidates for APCI. The presence of one or more heteroatoms (O, N, etc.) in an analyte increases the potential for the CI reaction. As such, many drugs, pesticides, herbicides, fungicides and natural products can be analyzed with APCI. Highly charged, thermally unstable and large molecular weight analytes should be avoided. These compounds are unacceptable as analytes and should also be avoided as matrix components since the lack of volatility and degradation products may increase maintenance of the APCI vaporization assembly and chamber.

Installation of APCI Chamber



1. Remove the Electrospray Ionization (ESI) chamber, if installed.

WARNING:
BURN HAZARD

Source parts become hot during operation. Use caution when installing or removing the API Source.


2. Connect the gas lines to the nebulizer assembly.
Nebulizing gas – top fitting.
Auxiliary gas - bottom fitting.


3. Connect the large (14-pin) round connector to the socket marked "API Housing."
4. Connect the small (4-pin) round connector to the socket marked "APCI Plug."

If this is the first time an APCI chamber is being installed on this instrument, do the following.

1. Change the instrument to APCI mode, by entering **APCI** on the command line.
2. Before connecting the corona high voltage connector, calibrate the needle current readback.
3. Select **Inst setup/others...** from the **Tools** menu.

4. Click Zero the APCI Needle Current, and click Execute.

5. After the routine is done, plug the corona high voltage connector into the socket marked "Needle Voltage."
6. Close and latch the chamber.

 WARNING: INHALATION HAZARD	The Source must be vented to an external fume hood or an external exhaust. Failure to properly vent the ion source can result in hazardous vapors being released into the laboratory.
--	---

APCI Daily Startup and Shutdown

Startup

1. Set the vaporizer and drying gas temperatures.
2. Click the spray icon  to turn on the gases and heaters.
3. When the vaporizer temperature has reached to a minimum of 75% of its set value, start the LC pump.
4. After waiting a few minutes for the temperatures to stabilize, the APCI is ready to use.

Shutdown

1. Turn off the LC pump flow. Allow the pump pressure to decrease and stabilize so flow into the interface stops.
2. Click the spray icon  to turn off the gases and heaters. The heaters switch off immediately, and after a few seconds, the gases switch off. The heaters will continue to cool.
3. Check that the mobile phase is no longer flowing before opening the spray chamber.

Tuning

Tune and mass axis calibrate of the 1200L mass spectrometer using the Electrospray Ionization (ESI) chamber. Refer to the Tuning section of the 1200L manual.

APCI Mobile Phase Selection

Most common HPLC mobile phases are ideal for LC/MS. Normal phase solvents such as dichloromethane, hexane and toluene, while unsuitable for ESI, are well suited to APCI. APCI can handle highly aqueous mobile phase compositions; however, you should use mobile phase additives with caution. Because APCI is a gas-phase ionization method and because mobile phase constituents are present in large excess compared to the sample, severe quenching of the analyte signal can occur. Particularly troublesome are additives such as formate and acetate in negative mode APCI. APCI usually uses a higher flow rate than ESI. HPLC flow rates for APCI should be between 0.1 to 2.0 mL/min, with optimum performance realized between 0.5 and 1.0 mL/min, supporting the most common 2.1 and 4.6 mm ID HPLC columns. When operating at higher or lower flow rates, the temperature of the APCI vaporizer should be raised or lowered appropriately.

Optimizing S/N for APCI

Do Optimization for each analysis. The detector response varies for each type of analyte and is dependant on the mobile phase composition and flow rate. Typical optimization normally involves flow injection or infusion of analyte into a mobile phase flowing at the analytical flow rate and composition while adjusting N₂ gas flows and temperatures to achieve maximum signal to noise. Use the syringe pump with flow rates less than 100 µL/min. However at very low flow rates there can be arcing to the shield.

The vaporizer temperature is the most significant parameter. Its optimum value is a function of:

- The sample - experimentally determined
- The mobile phase composition - more water requires a higher temperature, while a higher percentage of organic solvent requires a lower temperature
- The flow rate - higher flow rates require higher temperatures. If the temperature is very low, droplets will be observed at the outlet of the vaporizer tube. If the temperature is too high, the sample could thermally degrade.

Also important are the flow of the nebulizing and auxiliary gases, as well as the temperature and flow of the drying gas. Gas flows that are too low may cause an unstable baseline, while gas flows that are too high may reduce the size of the observed peaks.

The vaporizer can be moved back in the spray chamber to decrease the likelihood of fouling the vacuum capillary when the matrix contains a large amount of nonvolatile material (but at some expense in signal intensity). To adjust the position of the vaporizer relative to the capillary entrance, loosen the four mounting screws on top of the spray chamber and slide the vaporizer assembly forward or backward. Use the numbered scale on the right of the vaporizer assembly to assist in resetting the position. Tighten all four screws snugly when you are finished to avoid heat damaging the vaporizer. It must be firmly attached to the spray chamber for proper heat sinking.

Typical APCI Operating Parameters

WARNING:
BURN HAZARD

APCI Source may become hot during operation. Take appropriate precautions.

- **Drying Gas:** 12 psi (2 Liters/min) at 150 °C
- **Nebulizing Gas:** 60 psi (~1.5 liters/min)
Nitrogen in positive mode
Air in negative mode - Automatically switches if air is connected.
- **Auxiliary Gas:** 17 psi (2 liters/min), as hot as possible without decomposing sample (maximum temperature setting is 550 °C)
- **Corona Current:** 5-10 µA
- **Housing Temperature:** 50 °C

APCI Performance Evaluation

Operating Conditions for Testing the APCI

- **Drying Gas:** 12 psi at 150 °C
- **Nebulizing Gas:** 60 psi at 550 °C
- **Auxiliary Gas:** 17 psi
- **Corona Needle Current:** 5 μ A
- **Housing Temperature:** 50 °C
- **Solvent:** 50:50 methanol/water at 1 mL/min
- **Multiplier:** Set to optimum value
- **Scan Method:** 0.25 sec/scan, Mass peak width = 1.0, SIM width = 0.7
- **Reserpine:** Set the capillary between +80 to +120V and monitor m/z 609.3
- **4-nitrophenol:** Set the capillary between -45 to -65V and monitor m/z 138.0
- **Injection:** 5 μ L of either the reserpine or 4-nitrophenol solution, which is 2 pg/ μ L (10 pg injected.)

APCI Sensitivity Test

Positive APCI

Set the SIM scan mode to mass 609.3. A signal-to-noise of approximately 20:1 (RMS) or greater should be attained for 10 pg/ μ L reserpine (2 pg/ μ L, 5 μ L injection by FIA) at 1 mL/min. The mobile phase is 50:50 methanol/water. Check the results in Chro using 3 point smoothing. Right-click the peak to display the signal-to-noise ratio.

Negative APCI

Set the SIM scan mode to mass 138. A signal-to-noise of approximately 20:1 (RMS) or greater should be attained for 10 pg/ μ L nitrophenol (2 pg/ μ L, 5 μ L injection by FIA) at 1 mL/min. The mobile phase is 50:50 methanol/water. Check the results in Chro using 3 point smoothing. Right-click the peak to display the signal-to-noise ratio.

APCI General Maintenance

Cleaning

Clean the chamber daily with methanol to decrease the accumulation of contamination. If the chamber requires more aggressive cleaning, use slurry of aluminum oxide powder to polish the inside of the housing.

WARNING:
CHEMICAL HAZARD

Materials used for cleaning can produce hazardous vapors. Take appropriate precautions before beginning any cleaning procedure.

More extensive cleaning is required when sensitivity is lost or if excessive peak tailing occurs. Clean the drying gas spray shield, the vacuum capillary, the corona needle, and the quartz vaporizer tube.

NOTE: The quartz vaporizer tube may look dirty and the performance may not be affected. Replace it only if performance degrades.

Cleaning the Spray Shield

1. Remove the spray shield by gently rotating it off the mounting posts.
2. Wipe it with methanol or, if needed, scrubbing with an aluminum oxide powder/water slurry
3. Rinse with distilled water and then a methanol rinse.
4. Replace the spray shield by engaging one detent on the shield with one of the mounting posts and then using a gentle rotating motion to engage the other
5. Check that the capillary is visible and centered in the hole.

Cleaning the Vacuum Capillary


You can rinse the outer end of the capillary without venting the system.

Using a wash bottle spray a small amounts of water or methanol into capillary. You can remove the spray shield.

If more extensive cleaning of the inside of the capillary is required, the system must be vented, the capillary removed. Place it inside a graduated cylinder in an ultrasonic bath and clean it.

Cleaning the Corona Needle

Wipe the corona needle with methanol. If the deposits cannot be easily removed, replace the needle as follows.

Cutaway View of APCI Chamber

Replacing and Aligning the Corona Needle

Replace the corona needle needs if it is bent, damaged, or cannot be cleaned satisfactorily.

Remove the corona needle by grasping it with your fingers and pulling it straight out of its mount. **Be careful of the sharp point.** Insert a new or cleaned corona needle into the mounting hole and press firmly until it is fully seated.

The alignment should not change during normal needle cleaning or replacement. The needle is properly aligned when there is 1 cm between the end of the needle and the end of the vaporizer assembly.

While firmly holding the corona needle mount, loosen the set screw in the mount with a 1.5 mm hex key and rotate the needle mount until the needle is 1 cm below the vaporizer assembly, and then retighten the set screw. Do not remove the set screw from the mount. Do not remove the needle mount from the high voltage standoff. Hold the needle mount firmly while making adjustments.

Replacing the Nebulizing Needle

If the liquid needle becomes plugged, it can be easily replaced. The needle assembly is a length of 0.004 in. ID stainless steel tubing bonded to a length of 0.005 in. ID PEEK™ tubing.


Removal

Unscrew the fingertight PEEK fitting on top of the nebulizer assembly. Remove the old needle assembly from the top of the nebulizer assembly by pulling the PEEK tubing.

Installation

Carefully feed the new needle assembly (03-937544-03) into the hole in the nebulizer assembly, being careful not to bend the new needle. Slide a new PEEK fingertight ferrule and nut (28-211542-00) over the end of the PEEK tubing.

While ensuring that the PEEK tube is bottomed out, gently tighten the fitting.

Replacing the APCI Quartz Tube

It is possible for the quartz vaporizer tube to look dirty and not affect performance. Replace the quartz tube only if performance degrades.

Keep the components as clean as possible to avoid contaminating the interface.

 CAUTION	Be very careful while handling the quartz tube because it is very fragile.
 WARNING: BURN HAZARD	Components become very hot and can present a burn hazard. Ensure that all parts are cool before beginning this procedure.
 WARNING: CHEMICAL HAZARD	The spray chamber may contain hazardous chemicals. Use appropriate protection when handling.

Removing the Quartz Tube

1. Cool the APCI vaporizer.


NOTE: To cool the vaporizer, leave the spray icon on  and set the vaporizer temperature to zero. Leaving the gases on carries away residual heat.

2. Turn the gases and voltages off by clicking the spray icon  (blue LED on the front of the instrument must be off.)
3. Disconnect the gas lines from the bulkhead fittings on the front of the instrument.
4. Disconnect liquid connection to injector or column (not the fitting on top of the nebulizer.)
5. Disconnect the large round (14-pin) connector from the front panel of the instrument. The socket is marked "API Housing."
6. Disconnect the small (4-pin) round connector from the socket marked "APCI Heater."
7. With the spray chamber still closed and latched, loosen the two screws on top of nebulizer. (The nebulizer has the heat-sink fins.)
8. Remove nebulizer assembly (03-937326-01) by lifting straight up. Set it aside carefully so nebulizer tube does not get damaged or bent.
9. Carefully remove the black O-ring (03-930109-06) from the end of the quartz tube.

10. Remove the white ceramic insulating washer (03-937323-01) and set both aside.
11. Open the spray chamber. Leave it attached to the instrument by the mounting hinge.
12. Remove corona needle (03-937344-01) from its mount by pulling straight up. **Be careful of the sharp point.** Set aside.
13. After checking that the vaporizer has cooled, loosen and remove the threaded compression fitting (03-937340-01) from the end of the vaporizer tube inside the spray chamber.
14. Carefully pull the quartz tube out the bottom of the vaporizer (the end in the spray chamber). Be sure to keep the tube in line with the heater as you pull it out. Do not tilt it to the side. Use a slight twisting motion if the tube binds within the heater. The ferrule will remain attached to the tube. Discard both the quartz tube and the ferrule.

Installing the Quartz Tube

1. Close and latch the spray chamber.
2. Place the white ceramic insulating washer (03-937323-01) into the cutout on the vaporizer, covering the top of the heater. Gently press on the washer to seat it in the cutout.
3. Take the new quartz tube (03-937327-01) and place the O-ring (03-930109-06) over the end. The quartz tube should be barely protruding from the O-ring for this step.
4. Carefully lower the new quartz tube into the heating coil of the vaporizer. The O-ring on the quartz tube will be resting on the insulating washer when the quartz tube is properly installed. If the quartz tube does not go in all the way, GENTLY rock the insulating washer and/or the quartz tube with gentle pressure, taking care not to break the quartz tube. **DO NOT FORCE THE QUARTZ TUBE INTO THE HEATING COIL.** If the quartz tube snags the heating coil and damages it, the vaporizer assembly must be replaced. Once the quartz tube is in the heating coil, gently press on the end of the tube to push it slightly below the O-ring.

Proper Depth of Quartz Tube in O-ring

5. Open the spray chamber, keeping a finger on the quartz tube through Step 7 so it doesn't fall out.
6. Place a new polyimide ferrule (03-937339-01) over the exposed end of the quartz tube with the wider side of the ferrule toward the vaporizer assembly.
7. Replace the compression fitting and tighten finger tight. *NOTE: Excessive tightening can break the quartz tube.* The fitting is correctly tightened when gentle pressure on the end of the quartz tube cannot slide it within the vaporizer. Check to see that the O-ring is touching the insulating washer. The quartz tube is now fixed in place and no longer needs to be held by your finger.
8. Replace the corona needle by pressing it firmly into its mount until it bottoms.
9. Close and latch the spray chamber.
10. Replace the nebulizer assembly, lowering it straight down, using care to first pass the nebulizing tube into the quartz tube and then to engage the alignment pins (the gas lines should be on the right side). Excessive rocking of the nebulizer assembly can crack the quartz tube.
11. After checking to make sure the nebulizing assembly is straight, evenly tighten the two screws on top until snug.
12. Connect the large (14-pin) round connector to the socket marked "API Housing."
13. Connect the small (4-pin) round connector to the socket marked "APCI Heater."
14. Set the APCI temperature to 550 °C and click on the spray icon  to turn on the gases and heaters.
15. Allow the system to bake out for 15 minutes to remove surface contamination from the new quartz tube.
16. The APCI interface is ready to use. No adjustment is necessary if the corona needle mount has not been rotated. If the corona needle needs adjustment, follow the procedure above.

APCI Troubleshooting

If no peaks are observed, but the background spectrum looks normal, the temperature may be set too low for the sample.

If smaller than expected peaks are observed or the peak shape shows extreme tailing, the temperature may be set too low.

If the background spectrum has no ions, try increasing the corona current. If operating in negative mode, be sure the nebulizing gas is air.

If condensation is observed in the spray chamber, the housing temperature may be too low for the flow conditions or there has been insufficient warm-up time.

If the baseline is noisy, it may be as a result of poor nebulization. Check the gas flows and temperatures.

If the baseline is noisy and/or elevated, the chamber may need cleaning or the mobile phase may be contaminated.

If liquid drops are seen at the end of the quartz tube it may be because the mobile phase contains a high percentage of water flowing at a high flow rate. Increase the vaporizer temperature or reduce the flow rate until the drops are no longer visible.

If intense mass peaks are observed in the background spectrum, even without liquid flow, try cleaning, then baking out the spray chamber at maximum temperature for 2 hours. If the contamination still persists then the O-rings which seal the chamber and/or the chamber windows may be contaminated and should be replaced.

If the vaporizer temperature does not increase when turned on, contact your local Varian, Inc. representative for instrument repair and service.

Comparison of ESI and APCI Parameters

PARAMETER	ESI	APCI
Drying Gas	4L/min @ 350 °C	2L/min @ 150 °C
Nebulizing Gas	60 psi	60 psi
Auxiliary Gas	N/A	2L/min @ 550 °C
Inner Needle	Adjustable	Fixed
Needle Assembly Up/Down	Adjustable	Fixed
Needle Assembly Translation	Adjustable	Adjustable
Needle Voltage/Current	5000V	5 µA
Optimum Liquid Flow Rate	200-400 µL/min	0.5-1 mL/min
Housing Temperature	40 °C	50 °C

APCI Spare Parts

PART NUMBER	DESCRIPTION
03-937327-01	Quartz Vaporizer Tube
03-930109-06	Quartz Vaporizer Tube O-ring
03-937339-01	Quartz Vaporizer Tube Ferrule (Polyimide)
03-937344-01	Corona Needle
28-211542-00	PEEK 2-piece 10-32 Nut and Ferrule
03-937544-03	APCI Liquid Needle Assembly, 0.004 ID
03-937337-01	Spray Shield
27-402529-00	Spray Chamber O-ring
27-402690-00	O-ring, Spray Chamber Window
03-937323-01	Ceramic Insulating Washer (white)
03-920276-00	Aluminum Oxide, 600 Grit
88-999990-00	Applicator, Cotton Tipped, pkg. 100
03-937279-01	Reserpine Test Sample, 5 ng/μL, 6 Vials
03-937283-01	4-Nitrophenol Test Sample, 5 ng/μL, 6 Vials

These components only need replacing if they are damaged or lost.

Part Number	Description
03-937340-01	APCI Vaporizer Compression Fitting
03-937326-01	Nebulizer Assembly