

Agilent Nanospray/Nanodapter

Quick Start Guide

Installation 3

- Step 1. Set up the LC system 3
- Step 2. Install the Sensirion Liquid Flow Viewer software 7
- Step 3. Set up the LC/MS and install the Nanospray source 10

Flow Configuration 13

- Direct Flow with Nano Column 14
- Nano Flow with Trap Column 15
- Nano Flow with Trap Column and Divert Tee 16
- Nano Flow with Trap Column and Divert Valve 17

Nanodapter Checkout 18

- Step 1. Purge the pump and wash station 18
- Step 2. Prepare peptide standard 19
- Step 3. Start the LC with no column 21
- Step 4. Set up the LC/MS 21
- Step 5. Install the needle and set up of the system 22
- Step 6. Add IRM solution 25
- Step 7. Set up the worklist 27
- Step 8. Run the method and worklist 27
- Step 9. Analyze the data 28

Operation Guidelines 30

- Operation Guidelines 30
- To prevent stability issues 32
- To test the Nanodapter at different flow rates 33



Troubleshooting	35
If you have flow path blockages	35
If you have unstable flow	36
If you have unstable spray	37
If you have poor chromatography	38
If you have sample carryover	38
If you have poor sensitivity	39
Parts	40
PEEK-coated fused-silica capillary	40
G1988-64003 Nanodapter	41
Ultra-Low Dispersion Kit	43

Installation

Follow these steps to install and configure your Nanospray/Nanodapter System.

If you are installing just the Nanospray Source, go to “[Step 3. Set up the LC/MS and install the Nanospray source](#)” on page 10.

Step 1. Set up the LC system

The Nanodapter is supported on the 1290 Infinity II LC system (with G7167B Multisampler) and the 1290 Infinity LC system (with G4226A Autosampler).

The *Nanodapter/Nanospray Installation and Maintenance Guide* can be found on the *Nanospray Customer Information Disc*. The topics contained in that guide can also be found in the animated *6000 Series LC/MS Maintenance Guide* (revision C or higher).

Table 1 Nanodapter and LC System Setup

Step	Refer to
1 Properly install and configure the LC modules. Make sure all LC modules are in good operational condition and are free of contaminants.	LC module <i>Installation Guides</i> .

2 Installation

Step 1. Set up the LC system

Table 1 Nanodapter and LC System Setup







Step	Refer to
<p data-bbox="358 336 803 493">2 Install the Ultra-low Dispersion Needle Seat Assembly. Needle Seat Assembly for G7167B Multisampler (left) and G4226A Autosampler shown below.</p> <div data-bbox="358 539 775 739"><p data-bbox="358 687 519 713">(p/n 5067-5963)</p><p data-bbox="615 713 775 739">(p/n 5067-5189)</p></div>	<p data-bbox="821 336 1269 395"><i>Nanodapter/Nanospray Installation and Maintenance Guide</i></p> 
<p data-bbox="358 786 803 907">3 Install Nanodapter brackets. Brackets for G7167B Multisampler (left) and the G4226A Autosampler shown below.</p> <div data-bbox="358 956 815 1081"></div>	<p data-bbox="821 786 1269 845"><i>Nanodapter/Nanospray Installation and Maintenance Guide</i></p> 


Table 1 Nanodapter and LC System Setup

Step	Refer to
<p data-bbox="358 336 782 395">4 Mount the Nanodapter to the LC system.</p> <ul data-bbox="396 413 801 586" style="list-style-type: none"><li data-bbox="396 413 696 439">• Mount the Valve Drive.<li data-bbox="396 456 793 515">• Install the High-pressure Valve Head.<li data-bbox="396 532 801 586">• Mount the Nanodapter over the Valve Head.	<p data-bbox="825 336 1268 395"><i>Nanodapter/Nanospray Installation and Maintenance Guide</i></p> 
<p data-bbox="358 786 739 845">5 Connect the pump to the Nanodapter with flex tubing.</p>	<p data-bbox="825 786 1268 845"><i>Nanodapter/Nanospray Installation and Maintenance Guide</i></p> 

2 Installation

Step 1. Set up the LC system

Table 1 Nanodapter and LC System Setup

Step	Refer to
6 Connect the Injection Valve (Port 1) to the Flow Sensor of the Nanodapter with a 35-cm PEEK-coated fused-silica tubing.	<i>Nanodapter/Nanospray Installation and Maintenance Guide</i> 





Step 2. Install the Sensirion Liquid Flow Viewer software

Step 2. Install the Sensirion Liquid Flow Viewer software



The Sensirion Liquid Flow Viewer software is included on the *Customer Information Disc*.


- 1 Open the *Customer Information Disc*.
- 2 Open the **Liquid_Flow_Sensors > SCC1_USB_WindowsDriver** folder.
- 3 Double-click **CDM v2.10.00 WHQL Certified**.
- 4 Click **Extract** in the **FTDI CDM Drivers** screen.
- 5 If you are told that the requested operation requires elevation, click **OK**.
- 6 If you get a **Program Compatibility Assistant** message, click **Close**.
- 7 If you get the **FTDI CDM Drivers** screen again, click **Extract**.
- 8 If you get the **User Account Control** message, click **Yes**.




2 Installation

Step 2. Install the Sensirion Liquid Flow Viewer software


9 When the installation completes, click **Finish**.




10 From the Customer Information Disc, open the **USB-RS485_Sensor_Viewer** folder.



11 Double-click **Sensirion_USB_RS485_Sensor_Viewer_V2.32**.



12 When the installation completes, click **Finish**.



13 Close all opened program and reboot the computer.

14 Plug the Flow Sensor USB cable into the computer.

Step 2. Install the Sensirion Liquid Flow Viewer software

15 On your computer, open the folder **c:\Program Files (x86)\Sensirion AG\USB_RS485_SensorViewer**.

16 Double-click **UsbSensorViewerGUI**.


17 In the Product Selection dialog box:

- For **Sensor Product**, select **Liquid Flow Sensor**.
- For **COM Hardware**, select **RS485/USB Sensor Cable**.
- For **COM Port Settings**, change to the appropriate com port.

18 Click **OK**.

19 Click **Run**.

Name	Date modified	Type	Size
x86	11/11/2016 8:06 AM	File folder	
DeviceDriverRoles.dll	12/17/2013 8:28 AM	Application extens...	23 KB
FileIO.dll	9/6/2013 2:32 PM	Application extens...	24 KB
Gurock.SmartInspect.dll	12/18/2012 1:14 PM	Application extens...	120 KB
lowkit.dll	12/19/2012 1:24 PM	Application extens...	71 KB
Nini.dll	12/18/2012 1:13 PM	Application extens...	55 KB
Sensirion.Calc.dll	9/6/2013 2:32 PM	Application extens...	51 KB
Sensirion.Utilis.dll	9/6/2013 2:32 PM	Application extens...	136 KB
SensorCable.dll	2/27/2014 10:25 AM	Application extens...	330 KB
SHDLC.dll	2/27/2014 10:25 AM	Application extens...	540 KB
UsbSensorViewer.dll	2/28/2014 10:01 AM	Application extens...	189 KB
UsbSensorViewerGUI	2/28/2014 9:51 AM	Application	197 KB
ZedGraph.dll	12/19/2012 1:24 PM	Application extens...	300 KB




2 Installation

Step 3. Set up the LC/MS and install the Nanospray source

Step 3. Set up the LC/MS and install the Nanospray source

Table 2 LC/MS and Nanospray Set up

Step	Refer to
1 Properly install and configure the Q-TOF LC/MS. Make sure the LC/MS instrument meets Sensitivity checkout specs with an AJS ion source.	<i>Q-TOF LC/MS System Installation Guide</i>
2 Install the auxiliary air line	<i>G1995A Low Background Site Preparation Kit Installation Guide</i>




3 Install the flow restrictor.	<i>Nanospray/Nanodapter Installation and Maintenance Guide or G1995A Low Background Site Preparation Kit Installation Guide</i>
---------------------------------------	---




Step 3. Set up the LC/MS and install the Nanospray source

Table 2 LC/MS and Nanospray Set up

Step	Refer to
4	Install the camera and power cord <i>Nanodapter/Nanospray Installation and Maintenance Guide</i>



5	Set the position of the Counter Electrode and LED camera light. <i>Nanodapter/Nanospray Installation and Maintenance Guide</i>
---	--




2 Installation


Step 3. Set up the LC/MS and install the Nanospray source

Table 2 LC/MS and Nanospray Set up

Step	Refer to
6 Assemble the Infusion Sample Delivery System.	<i>Nanodapter/Nanospray Installation and Maintenance Guide</i>



- 7** Install the Needle Sleeve Assembly into the Nanospray Ion Source and adjust the needle position.
- Nanodapter/Nanospray Installation and Maintenance Guide*




Flow Configuration

Use the PEEK-coated fused-silica tubings and flex tubings to configure the flow of your Nanoflow system.

The Tee Adapters on the Nanodapter can be used to hold Trap Columns or a Tee Diverter in place.

Instructions to connect the correct tubings are shown in the animated *Nanospray/Nanodapter Installation and Maintenance Guide*, which can be found on the *Customer Information Disc*.



3 Flow Configuration


Direct Flow with Nano Column

Direct Flow with Nano Column

Use this configuration to check out the Nanodapter setup.

In this configuration, the sample flows directly from the Injection Valve to the Nano Column.

Nano Flow Direct



Nano Flow with Trap Column

In this flow configuration, the sample flows through a Trap Column before it goes into the Nano Column.


The length of the transfer PEEK-coated fused-silica capillary depends on the length of the Nano Column.

Recommended transfer capillary and column lengths are listed in [Table 3](#).

Table 3 Recommended capillary and column length combinations with Trap Column

Transfer PEEK-coated fused-silica tubing length	Nano Column length
35 cm	15 cm
20 cm	25 cm
10 cm	50 cm

Nano Flow with Trap Column



3 Flow Configuration

Nano Flow with Trap Column and Divert Tee

Nano Flow with Trap Column and Divert Tee

In this flow configuration, the sample flows through a Trap Column and then splits between waste (via the Divert Valve) and the Nano Column.


The length of the transfer PEEK-coated fused-silica capillary depends on the length of the Nano Column.

Recommended transfer capillary and Column lengths are listed in [Table 4](#)

Table 4 Recommended capillary and column length combinations with Trap Column

Transfer PEEK-coated fused-silica tubing length	Nano Column length
35 cm	15 cm
20 cm	25 cm
10 cm	50 cm

Nano Flow with Trap Column and Divert Tee



Nano Flow with Trap Column and Divert Valve

In this flow configuration, the sample flows through a Trap Column and then splits between waste (via the Divert Valve) and the Nano Column.


The length of the transfer PEEK-coated fused-silica capillary depends on the length of the Nano Column.

Recommended transfer capillary and Column lengths are listed in [Table 4](#)

Table 5 Recommended capillary and column length combinations with Trap Column

Transfer PEEK-coated fused-silica tubing length	Nano Column length
50 cm	15 cm
35 cm	25 cm
10 cm	50 cm

Nano Flow with Trap Column and Divert Valve



4 Nanodapter Checkout

Step 1. Purge the pump and wash station

Nanodapter Checkout

Step 1. Purge the pump and wash station

Before you begin, make sure you have:

- LCMS-grade acetonitrile
- LCMS-grade water
- 0.1% formic acid
- isopropyl alcohol

1 Prepare and install:

- Channel A: LCMS-grade water with 0.1% formic acid
- Channel B: 90:10 LCMS-grade acetonitrile:LCMS-grade water with 0.1% formic acid
- Wash solvent S1: 2.7% LCMS-grade acetonitrile:LCMS-grade water with 0.1% formic acid
- Seal wash solvent: 50:50 isopropyl alcohol:LCMS-grade water

2 Purge the pump and the wash solvent station.

Step 2. Prepare peptide standard

This step creates a 20 fmol/μL working solution that contains 20 fmol/μL HSA peptide standard mix in 0.1% formic acid, 15% LCMS-grade acetonitrile:LCMS-grade water.

Before you begin, make sure you have:

- LCMS-grade acetonitrile
- LCMS-grade water
- high-quality polypropylene centrifuge tube
- HSA peptide standard mix (p/n G2455-85001)
- 0.1% formic acid
- micropipette
- sample vial (p/n 9301-0977)
- vial cap (p/n 5182-3458)

- 1 Make a Dilution Solvent from the compounds in Table 6 in a high-quality polypropylene centrifuge tube.

Table 6 Dilution Solvent

Reagent	Percentage (volume)
LCMS-grade acetonitrile	15%
LCMS-grade water with 0.1% formic acid	85%

Do not use polyethylene labware because of the presence of extractable compounds. Do not use glassware to avoid formation of sodium adducts of the peptides.

- 2 Allow the vial of HSA peptide standard mix to come to room temperature.

4 Nanodapter Checkout

Step 2. Prepare peptide standard

- 3 Make a 1 pmol/ μL **Stock Solution**. Use a **micropipette** to measure the amounts listed in **Table 7**.
 - a Add 500 μL **Dilution Solvent** to the vial of **HSA peptide standard mix** to create the **Stock Solution**.
 - b Close the vial cap and mix the vial briefly on a vortex mixer.
 - c Allow the **Stock Solution** to stand for four to five minutes at room temperature to ensure dissolution of the **HSA peptide standard mix**.
 - d Mix briefly again on a vortex mixer to homogenize the **Stock Solution**. Do not overmix on the vortex mixer or the sample can become foamy.

Table 7 Stock Solution

Reagent	Amount (volume)
Dilution Solvent	500 μL
HSA peptide standard mix	vial

- 4 Make a 20 fmol/ μL concentration **Working Solution**. Do not use glass vials. Hydrophobic peptides will be lost on glass surfaces.
 - a Add 196 μL **Dilution Solvent** to a **sample vial**.
 - b Add 4 μL **Stock Solution** to the **Dilution Solvent** in the **sample vial**
 - c Close the **sample vial** with a **vial cap** and briefly mix on a vortex mixer.

Table 8 Working Solution

Reagent	Amount (volume)
Dilution Solvent	196 μL
Stock Solution	4 μL

Make sure no bubbles form in the **Working Solution**. Bubbles injected into the LC can cause variability in responses.

- d Immediately refrigerate the samples at 4°C, either in the refrigerator or in the Multisampler.

The samples must be analyzed within one day.

Step 3. Start the LC with no column

- 1 Open the **Sensirion Liquid Flow Viewer**. Run the program to display the split flow rate.
- 2 Set the pump flow to **0.06 mL/min** and turn on the pump.
With no column installed, the split flow stabilizes at 5000 nL/minute or greater.
- 3 Turn off the pump and allow the pump pressure to fall to less than 10 bar.

Step 4. Set up the LC/MS

- 1 Set the capillary voltage to 900 V. Leave the LC/MS in **Standby** mode.
- 2 Turn on the video capture program to view the needle position inside the Nanospray source.
- 3 Set the vertical setting on the Nanospray source to **12.0 mm** height.
This position puts the needle tip above the opening in the spray cap.
- 4 Set the horizontal setting on the Nanospray source to **3.0 mm** away from the front face of the spray cap.
This position puts the needle about halfway between the electrodes, but slightly closer to the spray cap on the left side of the picture.

4 Nanodapter Checkout

Step 5. Install the needle and set up of the system

Step 5. Install the needle and set up of the system

- 1 Insert the Needle Assembly (with needle and Nano column installed) into the Nanospray source.

Refer to the animated *Nanospray/Nanodapter Installation and Maintenance Guide* for instructions to assemble the Needle Assembly.

- 2 Turn on the pump. Set the flow rate to 0.083 mL/minute.

The split flow settles at around 300 nL/min. A drop of mobile phase forms in approximately one minute on the end of the needle. If neither of these events happen, locate and fix the leak before you proceed.

CAUTION

Use care when you type the pump flow rate. A value of "0.83" mL/minute instead of "0.083" mL/minute puts the system at very high pressure and flow. Flow rates greater than 400 nL/min can damage the needle by splitting the tip. This type of damage is difficult to see or detect. It can result in an erratic spray and the need for very high capillary voltages.

- 3 Adjust the pump flow rate up or down to achieve a split flow rate close to **300 nL/min**.

A suggested deviation is within 5%, +/- 15 nL/minute of the target value. Start at 0.050 mL/minute with a maximum flow rate of 0.150 mL/minute. Do not exceed 750 Bar back pressure set in the method.

CAUTION

To avoid over pressure and permanent damage to the Sensirion Flow Sensor, always set the LC pump pressure limit and flow response rate to:

- Pressure Limit: 750 Bar
- Flow Ramp up: 25 mL/minute
- Flow Ramp down: 25 mL/minute

CAUTION

Do not touch the column or capillary to any surface as column or capillary can become fractured.

Do not touch the needle tip to any surface as the needle tip can be damaged.

- 4 Insert the column and holder into the Nanospray source.

Step 5. Install the needle and set up of the system

The actual settings of the dials on the Nanospray source depends on the instrument and on how the back opening of the needle is trimmed during its installation. A typical view of the needle position is shown in [Figure 1](#).

The ideal needle position is a balance between robustness and sensitivity. You can achieve higher signal when the needle is close to the spray cap (left side) or in a low position (into the drying gas stream). But these positions shorten the life of the needle.



Figure 1 Typical needle position

5 Turn on the LC/MS.

A signal might not be visible if a drop is on the end of the needle.

6 Increase the capillary voltage in 50 V steps until you see a signal.

Typically, you need not apply more than 1300 V capillary voltage to see a signal.

4 Nanodapter Checkout

Step 5. Install the needle and set up of the system

7 Allow the signal to stabilize.

In general, higher capillary voltage gives greater signal but shortens the life of the needle. Capillary voltages over 1500 V are uncommon and difficult to control.

If the spray does not start:

- a** Turn off the pump and allow the drop to evaporate.
- b** As soon as the drop disappears, turn the pump on again.


8 Allow the system a few minutes to stabilize. Fine-tune the signal as needed:

- Adjust the Vernier positions by +/- 0.5 mm increments.
- Vary the capillary voltage also.
- After the needle is completely dried and the spray is established, a spray can be sustained with a lower capillary voltage. A lower capillary voltage can improve the life of the needle.

Step 6. Add IRM solution

- 1 Prepare the IRM Working Solution.


Refer to the animated *Nanospray/Nanodapter Installation and Maintenance Guide* for instructions.



- 2 Apply the IRM solution to the Nanospray source.

Refer to the animated *Nanospray/Nanodapter Installation and Maintenance Guide* for instructions.

Apply carefully. IRM compounds are difficult to remove if too much is applied.



- 3 You can adjust the IRM abundances to a certain degree by adjusting the drying gas temperatures, which affects the responses.

4 Nanodapter Checkout

Step 6. Add IRM solution

NOTE

High background in the system significantly diminishes the response of peptides. Identify and remove any source of high background to acquire data of high sensitivity and high quality. Sources of background include:

- the exhaust system in the lab to which the LC/MS is vented
- the lab air, the drain bottle and hose, the MS source, mobile phases used in the LC, and contamination from previous samples injected into the system
- An example of a system with low background is shown below (the 1221.9906 m/z response is the IRM signal).




Figure 2 Example of a system with low background. m/z 1221.9906 is the IRM signal.

Step 7. Set up the worklist

Set up a worklist to follow good laboratory practice:

- Run at least 2 or 3 blanks before running the samples to condition the column and stabilize the spray and mass spectrometer.
- Run blank after each batch of samples (for example, control samples) to reduce carry overs and increases the accuracy of the results.

Running a blank after each batch of samples also increases the life span of the column.

Step 8. Run the method and worklist

- 1 Start the worklist.

The method loads the sample at 300 nL/minute from the sample loop onto the column. This takes 8 minutes after the Multisampler switches to mainpass mode. No data is acquired during this time.

After the loading is complete, the Multisampler changes to bypass mode, the pump gradient starts, and the LC/MS acquisition starts.

The method parameter files **AcquisitionMethod_ReportNanodapter6550.xls** and **AcquisitionMethod_ReportNanodapter6530.xls** are provided on the Customer Information Disc.

Step 9. Analyze the data

- 1 Open the MassHunter Qualitative Analysis program.
- 2 After the run completes, load the data file into MassHunter Qualitative Analysis.
- 3 Extract the EIC of the following peptide m/z values. Use a +/- 20 ppm window.

Table 9 HSA Ions and Charge States

HSA Ions	Charge State
686.287019	+2
464.250360	+2
575.311146	+2
547.317433	+3
637.648743	+3
671.821021	+2
581.636215	+3

- 4 Integrate the data. Use the **Agile2** integrator.

Figure 3 shows typical results.

The abundance value you get can differ from this example due to the needle position, capillary voltage, tune parameters, preparation of the HSA synthetic peptides standard, and other factors.

Look for these results:

- Seven peaks due to the seven peptides are present.
- The FWHM (Full Width Half Maximum) value for the most intense peptide LVNEVTEFAK, 575.3111 m/z needs to be less than 0.25 minutes.
- The standard deviation of the retention time of this peptide needs to be less than 1%.

Note that FWHM and standard deviation values are guidelines for installation only.

Note that:

- The last peptide, HPYFYAPELLFFAK 581.6362 m/z , is hydrophobic. Its abundance depends on the age of the sample.
- Sensitivity is instrument-dependent, so no sensitivity metrics are provided.




Figure 3 EIC of the extracted ion for HSA

Operation Guidelines

Instructions to connect the correct tubings are shown in the animated *Nanospray/Nanodapter Installation and Maintenance Guide*, which can be found on the *Customer Information Disc*.

Operation Guidelines


- When attaching or replacing the spray needle, do not over-tighten the black Needle Nut. Always hand-tighten.

CAUTION



If you overtightened the black Needle Nut, you can break the needle. Damage to the needle inside the ferrule is difficult to diagnose.

To properly replace a needle:

- Insert a new spray needle (blunt end first) through a new black conductive ferrule until the needle extends approximately 1 to 2 mm. Make sure the spray needle can easily slide through the black conductive ferrule.
- Carefully use your fingernail or a flat object to push the spray needle until the end of the needle is flush with the conductive ferrule.
- Hand-tighten the Needle Nut until the needle is snug but not tight. Turn the Needle Nut an additional 120 degrees to seal the junction.
- Wait for small droplets to form on the end of the needle tip. For better visibility, shine a light onto the tip. The small droplet will glow when present.
- Insert the spray needle, Needle Nut, and column into the Needle Assembly.



- Spray needles, when first installed, require a low starting voltage. For a flow rate of 300 nL/min and a good source electrical ground, the voltage setting for optimum spray is usually between 1250 and 1300 V.
 - a Start with the Vcap at 1200 V and insert the Needle Assembly into the source housing.
 - b Watch the camera image while you slowly lower the needle midway into the electrode elements.
 - c If the droplet is observed, set the Vcap to 1300 volts. The electrospray should now start.
 - d If the droplet is observed, but the spray does not start, stop the flow and wait for the droplet to evaporate or gently shake the needle to remove the droplet to induce the start of spray. If the tip is damaged, replace the needle.
 - e If no droplet is observed, and no ions are detected on the MS system, the tip is blocked. Replace the tip.
- Never run 8 μm to 10 μm Nanospray tips at a flow rate over 500 nL/min. The tip can crack and yield erratic spray and unusable data. If the liquid wicks up the side of the tip, remove the column and tip from the instrument and inspect for damage.
- Make sure to attach the grounding wire from the Nanospray source to the instrument chassis. Resistance greater than 10 ohm can result in spray instability.



To prevent stability issues

- To avoid spray needle blockage from sample particles:
 - Clean the sample with a Solid Phase Extraction (SPE) column or cartridge before injection. Detergents and tissue debris from samples must be removed to prevent clogs in the needle or tubings.
 - Use a trapping column before the analytical column to trap all particulates and insoluble materials.
 - Completely digest protein before running the sample. Undigested protein can plug the spray needle or column.
 - For protein digest samples, keep the acetonitrile concentration under 60%. Do not raise the acetonitrile to 90% at the end of the run.

In some LC methods, the analytical column is flushed by holding the gradient at 90% acetonitrile for 10 minutes or more to elute the hydrophobic peptides or large precipitants from the column. These peptides can cause high back pressure or blockage, which requires a change of the trapping or analytical column.

- Shavings or particles from the conductive black Nanospray needle ferrule can cause blockage in the spray needle tip. To verify presence of particles in the spray needle, look for black particles at the needle tip under a microscope. If shavings are present, replace the black ferrule and spray needle.
- If white particles are observed under a microscope in the spray needle tip, replace the column.

CAUTION

Columns can degrade when running at high flow velocities (> 6 cm/seconds at high pressures).

To test the Nanodapter at different flow rates

- 1 Connect the orange PEEK-coated fused-silica capillary to the output of the Flow Sensor. Put the other end of the capillary into a small waste container.
- 2 Purge the UHPLC Pump with 50:50 acetonitrile:water for 5 minutes.
- 3 Set the pump Pressure limit to 750 bar.

CAUTION

To avoid over pressure and permanent damage to the Sensirion Flow Sensor, always set the LC pump pressure limit and flow response rate to:

- Pressure Limit: 750 Bar
- Flow Ramp up: 25 mL/minute
- Flow Ramp down: 25 mL/minute

- 4 Set the LC pump flow rate to 0.050 mL/minutes.
- 5 Change the acetonitrile:water solution to 3:97.
- 6 Open the Sensirion Flow Viewer and make sure it is visible on your computer display.
- 7 Make sure that the pressure listed in the LC Icon panel shows over 180 bar.
- 8 In the Sensirion Flow viewer, make sure the Nano flow is over 500 nL/minute.
- 9 When the Sensirion Flow viewer shows a steady plot across the screen, take a note of the LC pump pressure and the Nano flow rate.
- 10 Change the LC pump flow setting from 0.050 mL/minute to 0.100 mL/minute.
- 11 In the Sensirion Flow viewer, make sure that the Nano flow rate has changed from approximately 500 nL/minute to approximately 1000 nL/minute (or approximately doubled).
- 12 Increase the flow rate from 0.100 mL/minute to 0.150 mL/minute.
- 13 Check that the pressure has increased proportionately. The pressure can be as high as 700 bar or more.
- 14 Check that the Nano flow has increased by about the same ratio.

Table 10 shows example results of this test. Make sure you wait five minutes between readings.

5 Operation Guidelines

To test the Nanodapter at different flow rates

NOTE

Your actual readings can be different from those listed in [Table 10](#). To determine whether the Nanodapter system works properly, look for a proportional increase in LC pressure and Nano Flow rate when you increase the LC flow. When you increase the LC flow by two-fold, expect to see the LC pressure and Nano flow rate also increase by approximately two-fold.

Table 10 Example numbers to expect

UHPLC Pump Flow	UHPLC Pressure	Nano Flow Observed
0.050 mL/min	220 bar	550 nL/min
0.100 mL.min	439 bar	1,090 nL/min
0.150 mL/min	660 bar	1,620 nL/min
0.050 mL /min	219 bar	556 nL/min

Troubleshooting

If you have flow path blockages

Cause	Solution
1 Fused silica particles in system	<ul style="list-style-type: none"> • Solvent-rinse capillaries and fittings with isopropanol or acetonitrile prior to installation. • Avoid overtightening fittings. • Avoid excessive bending or coiling of the capillaries. If you coil to a radius of less than 40 mm, you will damage the capillaries. • Avoid kinking, bending or crushing capillaries with LC doors or cover panels. • Backflush capillaries if they are already blocked.
2 Particles from solvent or sample	<ul style="list-style-type: none"> • Always use LCMS-grade reagents. • Clean up samples per <i>Sample Preparation Guide</i>.
3 Nanospray needle blocked due to heat without flow	<ul style="list-style-type: none"> • If you do not plan to use the system for several days, remove the nanospray needle/column holder assembly and store in protective plastic sleeve. • If the needle is blocked, flush with 100% B1 for 30 min. This sometimes removes the blockage.
4 Nanospray needle or tubing blocked by detergent or sample tissue debris.	<ul style="list-style-type: none"> • Clean the sample before injection. Use the process appropriate to the nature of the sample and sample extraction method, such as SPE cartridges, filters or centrifugation.

5 Troubleshooting

If you have unstable flow

If you have unstable flow

Cause	Solution
1 Flow blockages (can be especially troublesome in direct injection mode with sampler in mainpass)	<ul style="list-style-type: none">• Replace tubing that is blocked.
2 System pressure too low	<ul style="list-style-type: none">• Keep system pressure higher than 20 bar at the pump outlet. If 20 bar cannot be achieved, leak-test the LC.

If you have unstable spray

Cause	Solution
1 Capillary voltage not set correctly	<ul style="list-style-type: none"> • Adjust the capillary voltage. <ul style="list-style-type: none"> • Using the 8 μm needle tips and liquid flow between 175 and 300 nL/minute, make sure the capillary voltage is between 1600 and 2000 volts. • New needles usually require less voltage, but need slightly more as the needle ages and the tip becomes eroded or enlarged.
2 Nanospray needle not positioned correctly	<ul style="list-style-type: none"> • See the animated <i>Nanospray/Nanodapter Installation and Maintenance Guide</i>.
3 Flow path blockages	<ul style="list-style-type: none"> • Replace the tubing that is blocked.
4 Nanospray needle tip is damaged or partly blocked. (You observe sputtering or split spray)	<ul style="list-style-type: none"> • Replace the needle, reinstall the needle/column holder assembly in the source, and adjust needle position.
5 Ferrule not making good seal with needle	<ul style="list-style-type: none"> • Replace the ferrule, reinstall the needle/column holder assembly in the source, and adjust needle position.
6 Flow too great or needle tip enlarged (You observe steady bowed stream of liquid.)	<ul style="list-style-type: none"> • Reduce the LC flow. • Replace the needle as above.

If you have poor chromatography

Cause	Solution
1 Gaps at LC connections	<ul style="list-style-type: none">• When you connect a capillary to a fitting or the column, push the capillary into the fitting firmly and smoothly to avoid gaps.• When connections are leaking, set column flow to zero, loosen the fitting, reinsert the fused silica and retighten the fitting. If you tighten the fitting without re-seating the fused silica tube, you may allow a gap to remain between the fused silica and the fitting. This will result in peak dispersion.
2 Note: Chromatography in enrichment column mode is generally not as good as in direct injection mode.	<ul style="list-style-type: none">• Switch to direct injection mode if that is an option. Consider:<ul style="list-style-type: none">• Your injection volume• Level of salts and other water-soluble contaminants.

If you have sample carryover

Cause	Solution
1 No needle wash	<ul style="list-style-type: none">• Set up needle wash for well-plate sampler.
2 Inappropriate needle wash solvent	<ul style="list-style-type: none">• Switch to a solvent combination in which your sample is completely soluble.
3 In enrichment column mode, residual hydrophobic peptides in the injection system	<ul style="list-style-type: none">• Inject several rounds of blanks with high organic until the peaks are gone.

If you have poor sensitivity

Cause	Solution
1 Detector gain adjustment needs to be redone	<ul style="list-style-type: none"> • See instructions for verifying detector setting in your LC/MSD Trap documentation.
2 Capillary voltage too high, causing corona which can destroy peptides	<ul style="list-style-type: none"> • Reduce the capillary voltage.
3 Sample degradation from sitting at room temperature	<ul style="list-style-type: none"> • Prepare fresh samples. • If you have the optional thermostat on the micro well-plate sampler, make sure it is turned on and set to 4°C.
4 Peptides adsorbed on vial surface	<ul style="list-style-type: none"> • Switch to a different vial material (glass or plastic).
5 Bad injection due to air bubble at bottom of vial	<ul style="list-style-type: none"> • Tap vial gently to dislodge air bubble.
6 Unstable spray	<ul style="list-style-type: none"> • See “If you have unstable spray” on page 37 of this chapter.

Parts

The Parts List for the Nanospray source are in the *Nanospray/Nanodapter Installation and Maintenance Guide* (animated) and the *LC/MS Maintenance Guide* (animated, revision C or later).

This section contains parts for the Nanodapter.

PEEK-coated fused-silica capillary



Each capillary comes with:

- 5043-0915 Mounting Tool
- 5-mm Seal-tight Screw ×2 (×1 in G1988-68063 and G1988-68068)
- 8-mm ID Front Ferrule ×2 (×1 in G1988-68063 and G1988-68068)

Table 11 PEEK-coated fused-silica capillary

Kit Part Number	Description
G1988-68063	3.8-m, beige
G1988-68068	3.2-m, orange
G1988-68069	50-cm, blue
G1988-68070	35-cm, blue
G1988-68071	70-cm, blue
G1988-68072	20-cm, blue
G1988-68073	10-cm, blue

G1988-64003 Nanodapter



6 Parts
G1988-64003 Nanodapter




Table 12 Nanodapter Parts List


Number	Part
❶	G1988-00003 Valve Mount Base
❷	G1988-68068 3.2-m orange PEEK-coated fused-silica capillary
❸	8121-2817 2-m SCC1-USB Sensor cable
❹	G1988-20051 Knob
❺	G1988-20060 Tee Adapter
❻	G1988-80016 Divert Tee
❼	5021-1823 400-mm Flex Tubing with 0100-2259 Nut and 0100-2258 Ferrule
❽	G1988-00006 Capillary Basket
❾	G1988-68063 3.8-m beige PEEK-coated fused-silica capillary
❿	G1988-60020 UHPLC Tee Clamp Assembly
①	0960-3228 Flow Pressure Sensor (underneath Capillary Basket)
②	5065-9925 15-cm 2.7 μm Nano Column
③	G1988-60006 Needle Sleeve Assembly
④	G1988-60028 Nanospray Source Tray Assembly
⑤	G1988-00007 Tray Cover

Ultra-Low Dispersion Kit

The Ultra Low Dispersion Kit is available for the G7167B Multisampler (left) and G4226A Autosampler.

Table 13 Ultra Low Dispersion Kits

Module	Kit Part Number	Seat Assembly Part Number
G7167B Multisampler	5067-5963	G4267-87020



G4226A Autosampler	5067-5189	G4226-87030
--------------------	-----------	-------------



www.agilent.com

In This Book

This guide contains information to get started with the Nanospray ion source and the Nanodapter accessory.

Detailed installation and maintenance instructions are found in the animated *Nanospray/Nanodapter Installation and Maintenance Guide* on the *Customer Information Disc*.

© Agilent Technologies, Inc. 2017

Revision A, January 2017



G1988-90001

Revision A0



Agilent Technologies