

BROOKFIELD KU-3

Viscometer

Operating Instructions

Manual No. M04-243



INSTRUMENTATION & SPECIALTY CONTROLS DIVISION

11 Commerce Boulevard, Middleboro, MA 02302 USA
Tel: 508-946-6200 or 800-628-8139 (USA excluding MA)
Fax: 508-946-6262 Internet <http://www.brookfieldengineering.com>

Table of Contents

| | |
|--|-----------|
| I. INTRODUCTION | 5 |
| I.1 Components | 5 |
| I.2 Options | 6 |
| 1.2.1 Optional Paste Spindle | 6 |
| I.3 Specifications | 6 |
| I.4 Sample Container Specifications | 8 |
| I.5 Utilities | 8 |
| I.6 Safety Symbols and Precautions | 8 |
| I.7 Viscosity Units | 8 |
| I.8 Cleaning | 9 |
| | |
| II. OPERATION | 10 |
| II.1 Set-up | 10 |
| II.2 Key Functions | 10 |
| II.3 Taking Measurements | 11 |
| II.4 Troubleshooting | 12 |
| | |
| Appendix A - KU-3 Calibration Information | 13 |
| Appendix B - Instrument Dimensions | 19 |
| Appendix C - References | 20 |
| Appendix D - Online Help and Additional Resources | 21 |
| Appendix E - Warranty Repair and Service | 22 |

I. INTRODUCTION

The Brookfield KU-3 Viscometer measures fluid viscosity in Krebs units. A paddle type spindle is driven at 200 RPM by a constant speed motor. The reaction torque of the spindle rotating at 200 RPM is converted to viscosity in Krebs Units. The digital display of the Viscometer shows viscosity in Krebs Units (KU), the associated grams value (g), and viscosity in centipoise (cP). The value for centipoise is a conversion from the Krebs value as described in the ASTM standard D562. The Viscometer will measure viscosity from 40 KU to 141 KU, at weights of 32 to 1,099 grams (the equivalent centipoise range is 27 - 5,274 cP). Application reference information can be found in ASTM D562.

I.I Components

The KU-3 Viscometer package includes:

- (1) KU-3 Viscometer, upright rod and base
- (1) Power Supply (AV-6)
- (1) Paddle spindle (KU-1030)
- (1) Power Cord, 115/230 V (P/N varies) (not shown)
- (1) Adapter for US pint and US 1/2 pint cans (KU-1004)
- (1) Operating instructions (M04-243)




Figure I-1

Please check to be sure that you have received all components, and that there is no damage. If you are missing any parts, please notify AMETEK Brookfield or your local authorized dealer immediately. Any shipping damage must be reported to the carrier.

Note: Keep all packing materials for future use.

I.2 Options

1.2.1 Optional Paste Spindle

The paste spindle, Part No. KU-1031 (Figure I-2), is a special spindle; it is not included in a standard shipment/order. The design consists of offset rod-type vanes, approximately 2 mm diameter x 19 mm long. This spindle is suitable for use with high consistency materials such as roller mill pastes. *Do not use this spindle to report normal Krebs Unit measurements or centipoise.* Record the grams value shown in the digital display and note that the paste spindle was used to make the viscosity measurement.




Figure I-2

I.3 Specifications

| | |
|----------------------------|---|
| Range: | 40–141 KU 32–1,099 g 27 - 5,274 cP* |
| Accuracy: | ± 1% of full scale range |
| Repeatability: | ± 0.5% of full scale range |
| Paddle speed: | 200 rpm ± 0.1 rpm |
| Net weight: | 18 lb. (8.2 kg) |
| Gross weight: | 22.5 lb. (10.2 kg) |
| Dimensions: | 21 x 12 x 17in - Carton |
| Operating Environment: | 0°C (32°F) to 40°C (104°F) 20% - 80% R.H.: non-condensing atmosphere |
| Ingress Protection Rating: | IP20 |

Electrical Certifications:

Conforms to CE Standards:

| | |
|------------------|--|
| BSEN 61326: | Electrical equipment for measurement, control and laboratory use - EMC requirements. |
| BSEN 61010-1: | Safety requirements for electrical equipment, for measurement, control and laboratory use. |
| BSEN 50581:2012: | Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances (ROHS). |

*Centipoise scale is for reference only. Do not use for calibration.

NOTICE TO CUSTOMERS:



This symbol indicates that this product is to be recycled at an appropriate collection center.

Users within the European Union:

Please contact your dealer or the local authorities in charge of waste management on how to dispose of this product properly. All AMETEK Brookfield offices and our network of representatives and dealers can be found on our website: www.brookfieldengineering.com.

Users outside of the European Union:
Please dispose of this product according to your local laws.





Figure I-3

The instrument can be identified by the item or model number. The item and model numbers appear on the Serial tag located at the back of the instrument. Below is an explanation on how to read these numbers.



- A. KU-3 Designation
- B. Standard Torque
- C. K: Magnetic Coupling
B: Ball Bearing Suspension
- D. Unused

I.4 Sample Container Specifications

Container Dimensions:



| | <u>Can Height</u> | <u>Bottom Flange Diameter</u> |
|----------------|-------------------|-------------------------------|
| U.S. Pint | 3.923" (9.97 cm) | 3.400" (8.64 cm) |
| U.S. Half Pint | 3.000" (7.62 cm) | 2.885" (7.33 cm) |
| U.S. Quart | 4.850" (12.32 cm) | 4.240" (10.77 cm) |

If you have a requirement for container dimensions other than those listed above, please contact AMETEK Brookfield or your local authorized dealer.

I.5 Utilities

VAC; Hz Limits: 100-240 VAC; 50/60 Hz $\pm 5\%$



Power Supply: 15 Watts, Class II Certified Plug-in Power Supply Rated: 12V@1.25VAC

-  Main supply voltage fluctuations are not to exceed $\pm 10\%$ of the nominal supply voltage.
-  Must be used with AV-6 power supply. Alternate power sources may cause damage to the instrument.





I.6 Safety Symbols and Precautions

Safety Symbols

The following explains safety symbols which may be found in this operating manual.

-  Indicates hazardous voltages may be present.
-  Refer to the manual for specific warning or caution information to avoid personal injury or damage to the instrument.

Precautions

-  If this instrument is used in a manner not specified by the manufacturer, the protection provided by the instrument may be impaired.
-  This instrument is not intended for use in a potentially hazardous environment.
-  In case of emergency, turn off the instrument and then disconnect the electrical cord from the wall outlet.
-  The user should ensure that the substances placed under test do not release poisonous, toxic or flammable gases at the temperatures to which they are subjected to during the testing.

I.7 Viscosity Units


The Krebs Unit is a special measure of viscosity that is not based on the Newtonian model of flow. The ASTM test method D562 was originally developed around the special conditions of an instrument that used gravity to drive a paddle spindle at 200 rpm. The weight required to achieve 200 rpm varied depending upon the viscosity of the fluid under test. The Krebs


unit was developed through the correlation of the weights used and the time required for 100 revolutions of the paddle. The ASTM standard also provided a correlation from Krebs Units to the scientific measure of viscosity in centipoise.

The Grams scale represents the weight required to drive the paddle through the test fluid at a rate of 200 rpm. The gravity drive system specified in ASTM D562 required that the weight be varied until 100 revolutions were achieved in 30 seconds (200 rpm). The KU-3 drives the unit at 200 rpm automatically and provides the grams value that would be required on the gravity drive system. The grams scale is not a viscosity unit.

The Centipoise scale is available through a correlation originally defined in the ASTM test method. Since this value is based on the Krebs Unit, it is not equivalent to centipoise values determined using other types of viscometers, such as the Brookfield DV1MRV. The centipoise values displayed on the KU-3 are for reference only. Comparisons to measured values from other instruments should not be made.

I.8 Cleaning


 Hands/fingers must be clean and free of residue sample. Not doing so may result in deposit build up on the upper part of the spindle shaft and cause interference between the shaft and the magnetic coupling.

 Be sure to remove the spindle from the instrument prior to cleaning. Severe instrument damage may result if the spindle is cleaned in place.

Instrument and Keypad: Clean with a dry, non-abrasive cloth. Do not use solvents or cleaners.

Note: Optional screen protector (Part No. KU-1404) is available from AMETEK Brookfield or your local authorized dealer.

Immersed Components (spindles): Spindles are made of stainless steel. Clean with a non-abrasive cloth and solvent appropriate for sample material.

 When cleaning, do not apply excessive force, which may result in bending spindles.

 Keep the recess in the viscometer base free from sample material.

II. OPERATION

II.1 Set-up

- 1) The KU-3 is supplied with the viscometer head mounted to the stand. Remove the entire viscometer from the packing container and place on a sturdy level surface.



Make sure the instrument is in a decent working environment (dust-free, moderate temperature, low humidity, etc.).



Make sure the instrument is on a level surface.

- 2) Move the operating handle to the top (upper most) position.
- 3) **Optional:** Attach the screen protector, Part No. KU-1404.
- 4) Plug the metal jacket of the power supply into the circular receptacle on the back of the viscometer. Screw the threaded collar on the metal jacket securely to the threads on the receptacle. Connect the power cord to the socket on the power supply and plug into the appropriate power source.
- 5) Insert the paddle spindle (KU-1030) into the viscometer spindle shaft. Align the cross bars of the spindle to the coupling. Insert the spindle shaft completely until the magnetic connection is achieved.
- 6) For **Quart Cans:** Place the quart can directly into the recess of the viscometer base.
- 7) For **Pint Cans:** Place the can adapter (KU-1004) into the recess of the viscometer base. The deep side of the adapter should be facing up. Place the pint can directly into the adapter.

Caution: Slide the pint can into the can adapter to avoid making contact with the paddle spindle. Contact between the can and paddle spindle may result in damage to the KU-3 Viscometer.

- 8) For $\frac{1}{2}$ **Pint Cans:** Place the can adapter (KU-1004) into the recess of the viscometer base. The deep side of the adapter should be facing down. Place the $\frac{1}{2}$ pint can directly into the adapter.

II.2 Key Functions

Figure II-1 shows the control keys on the face of the KU-3 Viscometer. The following describes the function of each key.



POWER KEY

Turns power On/Off to KU-3 Viscometer. Press the power key once to turn power ON. Press the power key twice to turn power OFF.



HOLD

Holds readings on the display. Press the Hold key during viscosity measurement to "Hold" the current reading. Readings will be presented on the display even when the handle is lifted and the motor stops.




Figure II-1




UP/DOWN ARROWS


Up and Down Arrow Keys return active to change display units. Change the display unit from KU → g → cP.

II.3 Taking Measurements

- 1) Pour sample into the container. Fill to 3/4 inch (20 mm) below the container lip.
- 2) Bring the sample in the container to the specified temperature.
- 3) Press the KU-3 power key (see Figure II.2).
- 4) Select the desired display (KU, g or cP) using the display arrows. (This setting may be changed at any time during the test.)
- 5) For *Quart Cans*: Place the sample container on the Viscometer base.
- 6) For *Pint & 1/2 Pint Cans*: Place the sample container directly on the can adapter mounted on the viscometer base.
- 7) Move the Viscometer handle down to the lowest position. This will automatically immerse the spindle into the fluid. If the correct amount of fluid has been put into the container, the surface of the fluid will be at the immersion mark of the spindle.

 **CAUTION:** When using the 1/2 pint can, do not lower the spindle directly into the container. The narrow diameter of the can requires the spindle to be introduced at an angle. Tilt the 1/2 pint can while lowering the viscometer.

- 8) The spindle will begin to rotate once the handle is within 1/2 inch of the lowest position.

 **CAUTION:** Do not touch the spindle or shaft while it is rotating. Touching the spindle during rotation may cause the KU-3 to over-range "EEEE". An over-range condition may damage the sensing system of the KU-3 viscometer

- 9) Wait 5 seconds for the display reading to stabilize. A display of "EEEE" indicates an over-range condition.
- 10) Press the Hold key. The display arrow keys may be used to change the display.
- 11) Raise the handle to the top position. This stops the spindle from rotating.

Note: When using 1/2 pint cans, you must lift the can and then tilt slightly to remove the spindle.

- 12) Pull down on the spindle to remove for cleaning.



Figure II-2

II.4 Troubleshooting

- 1) Display is frozen on a single reading:
 - Press the Hold key
 - Power Cycle KU-3

- 2) Handle movement is difficult:
 - Adjust the tension bolt (hex) on top of the mounting block (back of the KU-3)

Appendix A - KU-3 Calibration Information

The accuracy of the KU-3 is verified using viscosity standard fluids calibrated in Krebs Units, which are available from AMETEK Brookfield or your local authorized dealer. Note: Calibration should not be verified using the centipoise scale. The centipoise scale is for reference use only and is based on a correlation from the Krebs Unit details in the ASTM test method D562. Viscosity Standards are calibrated at 25°C. Available standards are listed in Table A-1 below:

| <u>Brookfield Viscosity Standard</u> | <u>Nominal Viscosity (KU)</u> | <u>Temperature (°C)</u> |
|--------------------------------------|-------------------------------|-------------------------|
| KU61 | 61 | 25 |
| KU73 | 73 | 25 |
| KU87 | 87 | 25 |
| KU99 | 99 | 25 |
| KU106 | 106 | 25 |

Brookfield Viscosity Standards

We recommend that Brookfield Viscosity Standard Fluids be replaced in accordance with the expiration date on the label. These fluids can be stored under normal laboratory conditions. Disposal should be in accordance with local, state and federal regulations. Material Safety Data Sheets are available on our website: www.brookfieldengineering.com/support/documentation.

Calibration Check Procedure

The frequency of the calibration check should be based on your company's standard practice for test and calibration of instruments.

- 1) Select any two viscosity standards listed in Table A-1 to perform your calibration check. The viscosity standards are very temperature sensitive. The viscosity value of the fluid will change with temperature so it is important to control the temperature to 25.0°C.

Note: Do not use viscosity standards calibrated in centipoise.

- 2) Pour the selected fluid into a standard US 1-pint can. The spindle and the fluid in the can should come to temperature equilibrium before proceeding with the calibration check.

It is important that the fluid and the spindle come to 25.0°C, ± 0.1°C, before proceeding with the calibration check.

- 3) Once the fluid and spindle have come to thermal equilibrium, place the sample container on the viscometer base, using the KU1-1004 adapter. Use the Arrow Keys to select "KU".
- 4) Move the viscometer handle down to the lowest position. The spindle will begin rotating once the handle is within 1/2 inch of the lowest position.
- 5) Wait five (5) seconds for the reading to stabilize. Press the Hold Key. **You will need to record measurements in KU and g.** These two values work together to let you interpret the calibration results.
- 6) After you have recorded the results in KU, press the Arrow Keys until "g" is displayed and record your reading in grams.

Interpretation of Calibration Check Results

When verifying the calibration of the KU-3, the instrument tolerance and viscosity standard fluid tolerance must be **combined** to calculate the total allowable error.

The KU-3 is accurate to ± 11 grams, which is 1% of the full scale range in grams. The Brookfield Viscosity Standard is accurate to $\pm 1\%$ of the stated viscosity in KU.

The total allowable error should be stated in KU. Since the instrument accuracy is stated in grams, you will have to use the comparison table (Table A-2) and convert from grams to KU.

Correct interpretation of your calibration results requires that you compare your readings in grams to the equivalent in KU. You must then bracket your reading with upper and lower limits based on the allowable error of ± 11.0 grams. Convert this acceptable range in grams (as defined by the upper and lower limits) to KU units.

Example: Calculate the allowable error of the KU-3 using fluid KU106; the stated viscosity of the fluid is 104.8 KU. The viscometer indicated a measured viscosity of 105 KU and 410 grams.

- 1) Measured results from the calibration check in grams were 410 grams. Locate 410 grams on the conversion chart (Table A-2).
- 2) The KU-3 is accurate to ± 11.0 grams. Starting from 410 grams, count 11 places above and below the 410 grams. This is called "bracketing" the acceptable range. In this case, the acceptable range will be from 399 grams to 421 grams.
- 3) Convert the acceptable range in grams to KU. Locate the minimum and maximum grams bracketed. Look to the right of each number for the conversion to KU. In this case, it will be 103.9 KU minimum and 105.7 KU maximum. The total difference between 103.9 and 105.7 KU is 1.8 KU. Therefore, the accuracy is ± 0.9 KU. ***This is the accuracy for the Viscometer in KU, ± 0.9 KU.***
- 4) Now that you have the accuracy for the instrument, you can add it to the accuracy of the fluid. The fluid is accurate to $\pm 1\%$ of the stated value in KU. The viscosity standard is calibrated at 104.8 KU, $\pm 1\%$ is equal to ± 1.0 KU.

$$\begin{array}{r} 0.9 \text{ KU (instrument accuracy)} \\ + \underline{1.0 \text{ KU (fluid accuracy)}} \\ \hline 1.9 \text{ KU (total allowable error)} \end{array}$$

- 5) Total allowable error for the calibration check in this example is 104.8 KU, ± 1.9 KU (102.9 KU to 106.7 KU). Since the measured reading of 105 KU falls within this range, the Viscometer is considered in calibration.




Table A-2 (continued)

| g | KU | cP | g | KU | cP | g | KU | cP | g | KU | cP | g | KU | cP | g | KU | cP |
|------------|--------------|-------------|------------|--------------|-------------|------------|--------------|-------------|------------|--------------|-------------|------------|--------------|-------------|------------|--------------|-------------|
| 626 | 121 | 2948 | 676 | 123.6 | 3194 | 726 | 126.6 | 3440 | 776 | 129.7 | 3686 | 826 | 132.6 | 3932 | 876 | 135 | 4177 |
| 627 | 121.1 | 2953 | 677 | 123.7 | 3199 | 727 | 126.6 | 3445 | 777 | 129.8 | 3691 | 827 | 132.6 | 3936 | 877 | 135 | 4182 |
| 628 | 121.1 | 2958 | 678 | 123.7 | 3204 | 728 | 126.7 | 3450 | 778 | 129.8 | 3695 | 828 | 132.7 | 3941 | 878 | 135.1 | 4187 |
| 629 | 121.2 | 2963 | 679 | 123.8 | 3209 | 729 | 126.7 | 3454 | 779 | 129.9 | 3700 | 829 | 132.7 | 3946 | 879 | 135.1 | 8192 |
| 630 | 121.2 | 2968 | 680 | 123.8 | 3213 | 730 | 126.8 | 3459 | 780 | 130 | 3705 | 830 | 132.8 | 3951 | 880 | 135.2 | 4197 |
| 631 | 121.3 | 2972 | 681 | 123.9 | 3218 | 731 | 126.9 | 3464 | 781 | 130 | 3710 | 831 | 132.8 | 3956 | 881 | 135.2 | 4202 |
| 632 | 121.3 | 2977 | 682 | 123.9 | 3223 | 732 | 126.9 | 3469 | 782 | 130.1 | 3715 | 832 | 132.9 | 3961 | 882 | 135.3 | 4207 |
| 633 | 121.4 | 2982 | 683 | 124 | 3228 | 733 | 127 | 3474 | 783 | 130.2 | 3720 | 833 | 132.9 | 3966 | 883 | 135.3 | 4212 |
| 634 | 121.4 | 2987 | 684 | 124 | 3233 | 734 | 127 | 3479 | 784 | 130.2 | 3725 | 834 | 133 | 3971 | 884 | 135.4 | 4217 |
| 635 | 121.5 | 2992 | 685 | 124.1 | 3238 | 735 | 127.1 | 3484 | 785 | 130.3 | 3730 | 835 | 133 | 3976 | 885 | 135.4 | 4222 |
| 636 | 121.6 | 2997 | 686 | 124.2 | 3243 | 736 | 127.2 | 3489 | 786 | 130.4 | 3735 | 836 | 133 | 3981 | 886 | 135.4 | 4227 |
| 637 | 121.6 | 3002 | 687 | 124.2 | 3248 | 737 | 127.2 | 3494 | 787 | 130.4 | 3740 | 837 | 133.1 | 3986 | 887 | 135.5 | 4232 |
| 638 | 121.7 | 3007 | 688 | 124.3 | 3253 | 738 | 127.3 | 3499 | 788 | 130.5 | 3745 | 838 | 133.1 | 3991 | 888 | 135.5 | 4236 |
| 639 | 121.7 | 3012 | 689 | 124.3 | 3258 | 739 | 127.3 | 3504 | 789 | 130.5 | 3750 | 839 | 133.2 | 3995 | 889 | 135.6 | 4241 |
| 640 | 121.8 | 3017 | 690 | 124.4 | 3263 | 740 | 127.4 | 3509 | 790 | 130.6 | 3754 | 840 | 133.2 | 4000 | 890 | 135.6 | 4246 |
| 641 | 121.8 | 3022 | 691 | 124.5 | 3268 | 741 | 127.5 | 3513 | 791 | 130.7 | 3759 | 841 | 133.3 | 4005 | 891 | 135.6 | 4251 |
| 642 | 121.9 | 3027 | 692 | 124.5 | 3272 | 742 | 127.5 | 3518 | 792 | 130.7 | 3764 | 842 | 133.3 | 4010 | 892 | 135.7 | 4256 |
| 643 | 121.9 | 3032 | 693 | 124.6 | 3277 | 743 | 127.6 | 3523 | 793 | 130.8 | 3769 | 843 | 133.4 | 4015 | 893 | 135.7 | 4261 |
| 644 | 122 | 3036 | 694 | 124.6 | 3282 | 744 | 127.6 | 3528 | 794 | 130.8 | 3774 | 844 | 133.4 | 4020 | 894 | 135.8 | 4266 |
| 645 | 122 | 3041 | 695 | 124.7 | 3287 | 745 | 127.7 | 3533 | 795 | 130.9 | 3779 | 845 | 133.5 | 4025 | 895 | 135.8 | 4271 |
| 646 | 122 | 3046 | 696 | 124.8 | 3292 | 746 | 127.8 | 3538 | 796 | 131 | 3784 | 846 | 133.6 | 4030 | 896 | 135.8 | 4276 |
| 647 | 122.1 | 3051 | 697 | 124.8 | 3297 | 747 | 127.8 | 3543 | 797 | 131 | 3789 | 847 | 133.6 | 4035 | 897 | 135.9 | 4281 |
| 648 | 122.1 | 3056 | 698 | 124.9 | 3302 | 748 | 127.9 | 3548 | 798 | 131.1 | 3794 | 848 | 133.7 | 4040 | 898 | 135.9 | 4286 |
| 649 | 122.2 | 3061 | 699 | 124.9 | 3307 | 749 | 128 | 3553 | 799 | 131.1 | 3799 | 849 | 133.7 | 4045 | 899 | 136 | 4291 |
| 650 | 122.2 | 3066 | 700 | 125 | 3312 | 750 | 128 | 3558 | 800 | 131.2 | 3804 | 850 | 133.8 | 4050 | 900 | 136 | 4295 |
| 651 | 122.3 | 3071 | 701 | 125.1 | 3317 | 751 | 128.1 | 3563 | 801 | 131.2 | 3809 | 851 | 133.8 | 4054 | 901 | 136 | 4300 |
| 652 | 122.3 | 3076 | 702 | 125.1 | 3322 | 752 | 128.2 | 3568 | 802 | 131.3 | 3813 | 852 | 133.9 | 4059 | 902 | 136.1 | 4305 |
| 653 | 122.4 | 3081 | 703 | 125.2 | 3327 | 753 | 128.2 | 3572 | 803 | 131.3 | 3818 | 853 | 133.9 | 4064 | 903 | 136.1 | 4310 |
| 654 | 122.4 | 3086 | 704 | 125.2 | 3332 | 754 | 128.3 | 3577 | 804 | 131.4 | 3823 | 854 | 134 | 4069 | 904 | 136.2 | 4315 |
| 655 | 122.5 | 3091 | 705 | 125.3 | 3336 | 755 | 128.4 | 3582 | 805 | 131.4 | 3828 | 855 | 134 | 4074 | 905 | 136.2 | 4320 |
| 656 | 122.6 | 3095 | 706 | 125.4 | 3341 | 756 | 128.4 | 3587 | 806 | 131.5 | 3833 | 856 | 134 | 4079 | 906 | 136.2 | 4325 |
| 657 | 122.6 | 3100 | 707 | 125.4 | 3346 | 757 | 128.5 | 3592 | 807 | 131.6 | 3838 | 857 | 134.1 | 4084 | 907 | 136.3 | 4330 |
| 658 | 122.7 | 3105 | 708 | 125.5 | 3351 | 758 | 128.6 | 3597 | 808 | 131.6 | 3843 | 858 | 134.1 | 4089 | 908 | 136.3 | 4335 |
| 659 | 122.7 | 3110 | 709 | 125.5 | 3356 | 759 | 128.6 | 3602 | 809 | 131.7 | 3848 | 859 | 134.2 | 4094 | 909 | 136.4 | 4340 |
| 660 | 122.8 | 3115 | 710 | 125.6 | 3361 | 760 | 128.7 | 3607 | 810 | 131.7 | 3853 | 860 | 134.2 | 4099 | 910 | 136.4 | 4345 |
| 661 | 122.8 | 3120 | 711 | 125.7 | 3366 | 761 | 128.7 | 3612 | 811 | 131.8 | 3858 | 861 | 134.3 | 4104 | 911 | 136.4 | 4350 |
| 662 | 122.9 | 3125 | 712 | 125.7 | 3371 | 762 | 128.8 | 3617 | 812 | 131.8 | 3863 | 862 | 134.3 | 4109 | 912 | 136.5 | 4354 |
| 663 | 122.9 | 3130 | 713 | 125.8 | 3376 | 763 | 128.9 | 3622 | 813 | 131.9 | 3868 | 863 | 134.4 | 4113 | 913 | 136.5 | 4359 |
| 664 | 123 | 3135 | 714 | 125.8 | 3381 | 764 | 128.9 | 3627 | 814 | 131.9 | 3872 | 864 | 134.4 | 4118 | 914 | 136.6 | 4364 |
| 665 | 123 | 3140 | 715 | 125.9 | 3386 | 765 | 129 | 3632 | 815 | 132 | 3877 | 865 | 134.5 | 4123 | 915 | 136.6 | 4369 |
| 666 | 123.1 | 3145 | 716 | 126 | 3391 | 766 | 129.1 | 3636 | 816 | 132 | 3882 | 866 | 134.5 | 4128 | 916 | 136.6 | 4374 |
| 667 | 123.1 | 3150 | 717 | 126 | 3395 | 767 | 129.1 | 3641 | 817 | 132.1 | 3887 | 867 | 134.6 | 4133 | 917 | 136.7 | 4379 |
| 668 | 123.2 | 3154 | 718 | 126.1 | 3400 | 768 | 129.2 | 3646 | 818 | 132.1 | 3892 | 868 | 134.6 | 4138 | 918 | 136.7 | 4384 |
| 669 | 123.2 | 3159 | 719 | 126.1 | 3405 | 769 | 129.3 | 3651 | 819 | 132.2 | 3897 | 869 | 134.7 | 4143 | 919 | 136.8 | 4389 |
| 670 | 123.3 | 3164 | 720 | 126.2 | 3410 | 770 | 129.3 | 3656 | 820 | 132.2 | 3902 | 870 | 134.7 | 4148 | 920 | 136.8 | 4394 |
| 671 | 123.3 | 3169 | 721 | 126.3 | 3415 | 771 | 129.4 | 3661 | 821 | 132.3 | 3907 | 871 | 134.8 | 4153 | 921 | 136.8 | 4399 |
| 672 | 123.4 | 3174 | 722 | 126.3 | 3420 | 772 | 129.5 | 3666 | 822 | 132.3 | 3912 | 872 | 134.8 | 4158 | 922 | 136.9 | 4404 |
| 673 | 123.4 | 3179 | 723 | 126.4 | 3425 | 773 | 129.5 | 3671 | 823 | 132.4 | 3917 | 873 | 134.9 | 4163 | 923 | 136.9 | 4409 |
| 674 | 123.5 | 3184 | 724 | 126.4 | 3430 | 774 | 129.6 | 3676 | 824 | 132.4 | 3922 | 874 | 134.9 | 4168 | 924 | 137 | 4413 |
| 675 | 123.6 | 3189 | 725 | 126.5 | 3435 | 775 | 129.7 | 3681 | 825 | 132.5 | 3927 | 875 | 134.9 | 4172 | 925 | 137 | 4418 |

Table A-2 (continued)

| g | KU | cP | g | KU | cP | g | KU | cP | g | KU | cP | g | KU | cP | g | KU | cP |
|------------|--------------|-------------|------------|--------------|-------------|-------------|--------------|-------------|-------------|--------------|-------------|-------------|--------------|-------------|-------------|------------|-------------|
| 926 | 137 | 4423 | 956 | 138.2 | 4571 | 986 | 139.4 | 4718 | 1016 | 140 | 4866 | 1046 | 140.5 | 5013 | 1076 | 141 | 5161 |
| 927 | 137.1 | 4428 | 957 | 138.3 | 4576 | 987 | 139.4 | 4723 | 1017 | 140 | 4871 | 1047 | 140.5 | 5018 | 1077 | 141 | 5166 |
| 928 | 137.1 | 4433 | 958 | 138.3 | 4581 | 988 | 139.4 | 4728 | 1018 | 140 | 4876 | 1048 | 140.6 | 5023 | 1078 | 141 | 5171 |
| 929 | 137.2 | 4438 | 959 | 138.3 | 4586 | 989 | 139.5 | 4733 | 1019 | 140.1 | 4881 | 1049 | 140.6 | 5028 | 1079 | 141 | 5176 |
| 930 | 137.2 | 4443 | 960 | 138.4 | 4591 | 990 | 139.5 | 4738 | 1020 | 140.1 | 4886 | 1050 | 140.6 | 5033 | 1080 | 141 | 5181 |
| 931 | 137.2 | 4448 | 961 | 138.4 | 4595 | 991 | 139.5 | 4743 | 1021 | 140.1 | 4891 | 1051 | 140.6 | 5038 | 1081 | 141 | 5186 |
| 932 | 137.3 | 4453 | 962 | 138.5 | 4600 | 992 | 139.6 | 4748 | 1022 | 140.1 | 4895 | 1052 | 140.6 | 5043 | 1082 | 141 | 5191 |
| 933 | 137.3 | 4458 | 963 | 138.5 | 4605 | 993 | 139.6 | 4753 | 1023 | 140.1 | 4900 | 1053 | 140.7 | 5048 | 1083 | 141 | 5195 |
| 934 | 437.4 | 4463 | 964 | 138.6 | 4610 | 994 | 139.6 | 4758 | 1024 | 140.1 | 4905 | 1054 | 140.7 | 5053 | 1084 | 141 | 5200 |
| 935 | 137.4 | 4468 | 965 | 138.6 | 4615 | 995 | 139.6 | 4763 | 1025 | 140.1 | 4910 | 1055 | 140.7 | 5058 | 1085 | 141 | 5205 |
| 936 | 137.4 | 4472 | 966 | 138.6 | 4620 | 996 | 139.7 | 4768 | 1026 | 140.2 | 4915 | 1056 | 140.7 | 5063 | 1086 | 141 | 5210 |
| 937 | 137.5 | 4477 | 967 | 138.7 | 4625 | 997 | 139.7 | 4772 | 1027 | 140.2 | 4920 | 1057 | 140.7 | 5068 | 1087 | 141 | 5215 |
| 938 | 137.5 | 4482 | 968 | 138.7 | 4630 | 998 | 139.7 | 4777 | 1028 | 140.2 | 4925 | 1058 | 140.7 | 5072 | 1088 | 141 | 5220 |
| 939 | 137.6 | 4487 | 969 | 138.8 | 4635 | 999 | 139.7 | 4782 | 1029 | 140.2 | 4930 | 1059 | 140.8 | 5077 | 1089 | 141 | 5225 |
| 940 | 137.6 | 4492 | 970 | 138.8 | 4640 | 1000 | 139.8 | 4787 | 1030 | 140.2 | 4935 | 1060 | 140.8 | 5082 | 1090 | 141 | 5230 |
| 941 | 137.6 | 4497 | 971 | 138.8 | 4645 | 1001 | 139.8 | 4792 | 1031 | 140.2 | 4940 | 1061 | 140.8 | 5087 | 1091 | 141 | 5235 |
| 942 | 137.7 | 4502 | 972 | 138.9 | 4650 | 1002 | 139.8 | 4797 | 1032 | 140.2 | 4945 | 1062 | 140.8 | 5092 | 1092 | 141 | 5240 |
| 943 | 137.7 | 4507 | 973 | 138.9 | 4654 | 1003 | 139.8 | 4802 | 1033 | 140.3 | 4950 | 1063 | 140.8 | 5097 | 1093 | 141 | 5245 |
| 944 | 137.8 | 4512 | 974 | 139 | 4659 | 1004 | 139.8 | 4807 | 1034 | 140.3 | 4954 | 1064 | 140.8 | 5102 | 1094 | 141 | 5250 |
| 945 | 137.8 | 4517 | 975 | 139 | 4664 | 1005 | 139.9 | 4812 | 1035 | 140.3 | 4959 | 1065 | 140.8 | 5107 | 1095 | 141 | 5254 |
| 946 | 137.8 | 4522 | 976 | 139 | 4669 | 1006 | 139.9 | 4817 | 1036 | 140.3 | 4964 | 1066 | 140.9 | 5112 | 1096 | 141 | 5259 |
| 947 | 137.9 | 4527 | 977 | 139.1 | 4674 | 1007 | 139.9 | 4822 | 1037 | 140.3 | 4969 | 1067 | 140.9 | 5117 | 1097 | 140 | 5264 |
| 948 | 137.9 | 4532 | 978 | 139.1 | 4679 | 1008 | 139.9 | 4827 | 1038 | 140.4 | 4974 | 1068 | 140.9 | 5122 | 1098 | 141 | 5269 |
| 949 | 138 | 4536 | 979 | 139.1 | 4684 | 1009 | 139.9 | 4832 | 1039 | 140.4 | 4979 | 1069 | 140.9 | 5127 | 1099 | 141 | 5274 |
| 950 | 138 | 4541 | 980 | 139.2 | 4689 | 1010 | 139.9 | 4836 | 1040 | 140.4 | 4984 | 1070 | 140.9 | 5132 | | | |
| 951 | 138 | 4546 | 981 | 139.2 | 4694 | 1011 | 139.9 | 4841 | 1041 | 140.4 | 4989 | 1071 | 140.9 | 5136 | | | |
| 952 | 138.1 | 4551 | 982 | 139.2 | 4699 | 1012 | 140 | 4846 | 1042 | 140.4 | 4994 | 1072 | 140.9 | 5141 | | | |
| 953 | 138.1 | 4556 | 983 | 139.3 | 4704 | 1013 | 140 | 4851 | 1043 | 140.5 | 4999 | 1073 | 140.9 | 5146 | | | |
| 954 | 138.2 | 4561 | 984 | 139.3 | 4709 | 1014 | 140 | 4856 | 1044 | 140.5 | 5004 | 1074 | 140.9 | 5151 | | | |
| 955 | 138.2 | 4566 | 985 | 139.3 | 4713 | 1015 | 140.0 | 4861 | 1045 | 140.5 | 5009 | 1075 | 140.9 | 5156 | | | |

Appendix B - Instrument Dimensions



Appendix C - References

The KU-3 Viscometer is compatible with:

ASTM D562 Standard Test Method for Consistency of Paint Using the Stormer Viscometer

Appendix D - Online Help and Additional Resources

www.brookfieldengineering.com**

The Brookfield website is a good resource for additional and self-help whenever you need it. Our website offers a selection of “how-to” videos, application notes, conversion tables, instructional manuals, material safety data sheets, calibration templates and other technical resources.

<http://www.youtube.com/user/BrookfieldEng>

Brookfield has its own YouTube channel. Videos posted to our website can be found here as well as other “home-made” videos made by our own technical sales group.

Viscosityjournal.com

Brookfield is involved with a satellite website that should be your first stop in viscosity research. This site serves as a library of interviews with experts in the viscosity field as well as Brookfield technical articles and conversion charts. Registration is required, so that you can be notified of upcoming interviews and events, however, this information will not be shared with other vendors, institutions, etc..

Article Reprints

- **Available in Print Only**
- Brookfield has an extensive library of published articles relating to viscosity, texture and powder testing. Due to copyright restrictions, these articles cannot be emailed. Please request your hardcopy of articles by calling our customer service department directly or by emailing: MA-MID.sales@ametek.com
- **Available Online**
- Brookfield has a growing number of published articles that can be downloaded directly from the Brookfield website. These articles can be found on our main site by following this path: [http://www.brookfieldengineering.com/support/documentation/article reprints](http://www.brookfieldengineering.com/support/documentation/article%20reprints)

More Solutions to Sticky Problems

Learn more about viscosity and rheology with our most popular publication. This informative booklet will provide you with measurement techniques, advice and much more. It’s a must-have for any Brookfield Viscometer or Rheometer operator. More Solutions is available in print and also as a downloadable pdf on the Brookfield website by following this path:

<http://www.brookfieldengineering.com/support/documentation>

Training/Courses

Whether it is instrument-specific courses, training to help you better prepare for auditing concerns, or just a better understanding of your methods, who better to learn from than the worldwide leaders of viscosity measuring equipment? Visit our Services section on our website to learn more about training.

** Downloads will require you to register your name, company and email address. We respect your privacy and will not share this information outside of Brookfield.

Appendix E - Warranty Repair and Service

Brookfield Viscometers are guaranteed for one year from date of purchase against defects in materials and workmanship. They are certified against primary viscosity standards traceable to the National Institute of Standards and Technology (NIST). The Viscometer must be returned to **AMETEK Brookfield** or the authorized dealer from whom it was purchased for a warranty evaluation. Transportation is at the purchaser's expense. The Viscometer should be shipped in the packaging originally provided with the instrument. If returning to Brookfield, please contact us for a return authorization number prior to shipping. Failure to do so will result in a longer repair time.

*For a copy of the Repair Return Form, go to the Brookfield website,
www.brookfieldengineering.com*

For repair or service in the **United States** return to:

AMETEK Brookfield
11 Commerce Boulevard
Middleboro, MA 02346 U.S.A.
Telephone: (508) 946-6200 Fax: (508) 923-5009
www.brookfieldengineering.com

For repair or service outside the United States, consult AMETEK Brookfield or the authorized dealer from whom you purchased the instrument.

For repair or service in the **United Kingdom** return to:

AMETEK (GB) Limited
Brookfield Technical Centre
Stadium Way
Harlow, Essex CM19 5GX, England
Telephone: (44) 1279/451774 Fax: (44) 1279/451775
www.brookfield.co.uk

For repair or service in **Germany** return to:

AMETEK GmbH
Hauptstrasse 18
D-73547 Lorch, Germany
Telephone: (49) 7172/927100 Fax: (49) 7172/927105
www.brookfield-gmbh.de

For repair or service in **China** return to:

AMETEK Commercial Enterprise (Shanghai) Co., Ltd Guangzhou Branch
Room 810 Dongbao Plaza, No. 767 East Dongfeng Road
Guangzhou, 510600 China
Telephone: (86) 20/3760-0548 Fax: (86) 20/3760-0548
www.brookfield.com.cn

On-site service at your facility is also available from Brookfield. Please contact our Service Department in the United States, United Kingdom, Germany or China for details.

