

High Performance Liquid Chromatograph

i-Series

Specifications



Maximum Reliability and Stability

— Fundamental functions assure analysis results —



Advanced



High Performance Liquid Chromatograph

Use of Multiple Detectors Expands Application Range

In addition to the UV-visible (UV/VIS) absorbance detector or photodiode array (PDA) detector included as standard, a fluorescence detector or differential refractive index detector can be added.

Excellent Baseline Stability Unaffected by Circumstances

The UV/VIS detector and the PDA detector employ dual-temperature control (TC-Optics and flow cell) and provide measurements with a stable baseline hardly affected by room temperature fluctuation.

Supports High-Speed Multi-Analyte Processing

A 14-second injection cycle maximizes the number of samples that can be processed.

Moreover, a total of 1536 samples can be accommodated in right and left sample racks.

Autosampler Enhances Data Reliability

Excellent reproducibility for injection volumes less than 1 μ L, wide linearity range and ultra-low carryover (<0.0025%) improve the reliability of data, especially for analyses of precious biological samples and direct analyses of concentrated samples.

Open Access Sample Placement

A direct access mechanism on sample racks allows the user to place the sample on racks that are not involved in sample injection even during analysis.

Furthermore, racks can be shared by multiple analysts, without interrupting the analysis of samples placed by others. Overall, this function enhances work efficiency.



System Monitoring via Smart Devices

System status and chromatograms can be viewed remotely from smart devices and home computers. Situations can be checked without needing to visit the laboratory.



Window on a Smartphone



Refined Usability

Control panel with a color LCD touch panel allows anyone to operate the instrument, regardless of experience level. Easily and reliably perform routine maintenance following on-screen instructions.

Displays Chromatogram in Real Time

The chromatogram real-time monitor allows the user to immediately confirm the success or failure of data, even in a computer-less laboratory environment.

Large-Capacity Column Oven Supports Up to 90°C

The Forced-air circulation method enhances column temperature stability. Maximum operating temperature of 90°C allows high-temperature analyses such as sugar analysis. Moreover, three 300 mm long columns or six 100 mm long columns can be accommodated.

Quaternary Solvent Delivery Unit

A 10 μ L micro plunger ensures accurate quaternary gradient delivery. Optional reservoir switching valve further extends the solvent selection to seven so that the solvent for the flow path rinsing can be set.

Auto Shutdown Function Reduces Power Consumption

After analysis is complete, the auto shutdown function minimizes power consumption in standby mode and can reduce power consumption by at least 95% compared to normal standby mode.

Compact Footprint

The i-Series brings together all the functions required for LC analysis in an integrated form. With its space-saving design, which is only 410 mm wide, three units can be installed on a laboratory bench compared to only two comparable instruments from other companies. The instrument footprint does not change even if another detector is installed.

Specifications

	Model	LC-2050 UV model without sample cooler	LC-2050C (UV model)	LC-2050C 3D (PDA model)	LC-2050C LT (detector-less model	LC-2060C (UV model)	LC-2060C 3D (PDA model)	LC-2060C 3D MT (PDA model)			
	P/N	S228-65820-58	5228-65820-58 \$228-65821-58 \$228-65822-58 \$228-65823-58 \$228-65824-58 \$228-65825-58 \$228-65826-								
	Degassing unit	Five Lines: Mobile phase 4 + Rinse solution 1 (Volume 400 μL)									
	Pumping method	Parallel-type double plunger									
	Pulsation			≤ 0.1 MPa	(1.0 mL/min,10 M	IPa, Water)					
	Flow rate setting range			0	.0001 to 10 mL/m	in					
	Flow rate accuracy	(0	or $\leq \pm 2 \mu L/min$, b 0.01 to 2 mL/min, $\leq \pm$ (2 to 5 mL/min, S ₁	Specified conditio :2%	$\leq \pm 1\%$ or $\leq \pm 2$ µL/min, below whichever is greater (0.01 to 2 mL/min,1 to 40 MPa, Specified condition) $\leq \pm 2\%$ or $\leq \pm 2$ µL/min, below whichever is greater (0.01 to 3 mL/min, 40 to 60 MPa, Specified condition)						
	Flow rate precision		<	0.06 %RSD or < 0	.02 minSD, below	whichever is grea	ter				
	Configuration			Four-sol	vent low-pressure	gradient					
Pump	Gradient / range of set concentrations			0 to	100%, in 0.1% s	teps					
	Gradient / concentration accuracy		±().5% (0.1 to 2 mL/	min, 1 to 20 MPa,	, Specified condition	on)				
	Gradient / concentration precision			±0.1% (1 mL/	min,10 MPa, Speci	fied condition)					
	Maximum pressure		50 MPa (0.000 22 MPa (5.000	1 to 5 mL/min) 1 to 10 mL/min)	44 MF	Pa (0.0001 to 3 mL Pa (3.0001 to 5 mL Pa (5.0001 to 10 m	/min)				
	System Delay Volume) μL ι μL, 1100 μL)		Ο μL Ι μL, 1100 μL)	505 μL (Flow path 1) 765 μL (Flow path 2)				
	Injection method	Total-volume sample injection									
	Injection volume accuracy	±1% (50 μL, N = 10)									
	Injection volume setting range	0.1 to 100 μL 0.1 to 50 μL (Option: 0.1 to 50 μL, 1 to 500 μL, 1 to 500 μL, 1						L, 1 to 2,000 μL)			
utosampler	Injection volume reproducibility	RSD < 0.20% (5.0–2000 μL) RSD < 0.25% (2.0–4.9 μL) RSD < 0.5% (1.0–1.9 μL) RSD < 1.0% (0.5–0.9 μL)									
sam	Cross-contamination	0.0025% (Caffeine, Specified condition)									
Auto	Injection cycle time	Min. 14 sec (Specified condition)									
_	Samples for processing		ı	336 (1 mL), 216	(1.5 mL), 112 (4 m	L), 4 (MTP/DWP)					
	Sample cooler	Not included			ture needs to be 3 m temperature nee		higher to go up to				
	Injection linearity			9999 ecified condition)	S (1 to 100 μL,	> 0.9999 µL, Standard Samp pecified condition > 0.9999 , 100 µL Sample Lo Specified condition	oop (Option),				
	Heating and cooling method			Force	d air circulation m	ethod					
_	Containable column size		6 pied	ces at 10 cm max.,	30 cm		2 pieces at 10 to 30 cm				
Column Oven	Temperature control range	Room temperature – 12 to 90°C, Setting range 4 to 90°C									
Colui	Temperature control precision				±0.1°C						
	Temperature stability	±0.8°C (Specified condition)									
	Flow rate switching valve			Max	. 1 рс			Not allowed			

	Model	LC-2050 (UV model (without sample cooler)	LC-2050C (UV model)	LC-2060C (UV model)				
	P/N	\$228-65820-58	\$228-65821-58	S228-65824-58				
	Wavelength range		190 to 700 nm					
	Spectral bandwidth	8 nm						
	Wavelength accuracy	≤ ±1 nm						
	Wavelength reproducibility	≤ ±0.1 nm						
	Noise level	≤ ±2.5 × 10 ⁻⁶ AU, (250 nm, Specified condition)						
tor	Drift	≤ 100 × 10 ⁻⁶ AU/h (250 nm, Specified condition)						
UV Detector	Simultaneous monitoring of 2 wavelengths	Enable (Any 2 wavelengths of 190 to 370 nm or 371 to 700 nm)						
	Linearity	Up to 2.5 AU (5%)						
	Sampling rate	Up to 100 Hz						
	Light source	Deuterium (D ₂) lamp						
	Flow cell	12 μL (10 mm	, TC), 12 MPa	8 μL (10 mm, TC), 12 MPa				
	Option cell		μL (10 mm, TC), 5 μL (5mm, TC)	Conventional: 12 μL (10 mm, TC), Semi-micro: 2.5 μL (5 mm, TC)				

	Model	LC-2050C 3D (PDA model)	LC-2060C 3D (PDA model)	LC-2060C 3D MT (PDA model)				
	P/N	S228-65822-58	\$228-65825-58	S228-65826-58				
	Wavelength range		190 to 800 nm					
	Spectral resolution	1.4 nm (Specified condition)						
	Slit width	1.2 nm, 8 nm						
	Device resolution	0.6 nm/pixel						
or	Number of photodiode array elements	1024						
Detector	Wavelength accuracy	≤±1 nm						
	Noise level	≤ ±3 × 10 ⁻⁶ AU (250 nm, reference: 350 nm, Specified condition)						
PDA	Drift	≤ 500 × 10 ⁻⁶ AU/h (250 nm, reference: 350 nm, Specified condition)						
	Linearity	Up to 2 AU (5%)						
	Sampling rate	Up to 100 Hz						
	Light source	Deuterium (D ₂) lamp (Standard), tungsten (W) lamp (option)						
	Flow cell	10 μL (10mm, TC), 12 MPa	8 µL (10 mm,	TC), 12 MPa				
	Option cell	High-Speed: 8 μL (10 mm, TC), Semi-micro: 2.5 μL (5 mm, TC)	Conventional: 10 μL (10 mm, TC), Semi-micro: 2.5 μL (5 mm, TC)					

Model		LC-2050 (UV model without sample cooler)	LC-2050C (UV model)	LC-2050C 3D (PDA model)	LC-2050C LT (detector-less model)	LC-2060C (UV model)	LC-2060C 3D (PDA model)	LC-2060C 3D MT (PDA model)	
P/N		S228-65820-58	S228-65821-58	S228-65822-58	S228-65823-58	S228-65824-58	S228-65825-58	S228-65826-58	
	Dimensions		W410×H605×D500 mm (Not including reservoir tray height)						
ST	Weight	58 kg	63 kg		53 kg	63 kg		64 kg	
eous	Available pH range	1 to 13							
Miscellar	Materials for parts in contact with liquids	Stainless steel (SUS316L, SUS316), FEP, PEEK, PTFE, perfluoroelastomer, ruby, sapphire, Hastelloy® C, GFP, ceramic, PFA, quartz, PPS Stainless steel (SUS316L, SUS316), FEP, PEEK, PTF perfluoroelastomer, ruby, sapphire, Hastelloy C, UHMWPE, ceramic, PFA, quartz, PPS							
	Workstation	LabSolutions™ LC/GC Ver.5.103 or later, LabSolutions DB/CS Ver.6.103 or later (Incompatible with LCsolution™)							

Optional Detector Specifications

RID-20A

	RID-20A (\$228-65306-58)
Reflective index measurement range	1 to 1.75 RIU
Noise level	≤ 2.5 nRIU
Drift	≤ 0.1 μRIU/h
Range	A mode: 0.01 to 500 μRIU P and L modes: 1 to 5000 μRIU
Response	No filtering, 0.05 to 10 sec, 11 steps
Polarity switching	With a switch
Zero adjustment	Auto zero, auto optical zero, baseline shift functions
Maximum operating flow rate	20 mL/min (150 mL/min with an option)
Temperature control of cell unit	30 to 60°C (0.01°C steps)
Cell capacity	9 μL
Material in contact with liquid	SUS316L, quartz, PTFE, Al ₂ O ₃ , ETFE
Maximum operating pressure	0.4 MPa (4 kgf/cm²)
Operating temperature range	4 to 35°C
Dimensions and weight	W260 × D420 × H140 mm, 12 kg

Note: Hexafluoroisopropanol (HFIP) cannot be used as the mobile phase.

RF-20A/RF-20Axs



	RF-20A (S228-65304-58)	RF-20Axs (S228-65305-58)			
Light source	Xenon lamp	Xenon lamp, low-pressure mercury lamp (To check wavelength accuracy)			
Wavelength range	0, 200 to 650 nm	0, 200 to 750 nm			
Spectral bandwidth	20	nm			
Wavelength accuracy	±2	nm			
Wavelength precision	±0.2	2 nm			
S/N	Water Raman peak S/N 1200 min. Water Raman peak S/N 2000 Low background S/N > 9000 Low background S/N > 120				
Cell capacity	SUS316L, PTFE (fluororesin), quartz				
Cell temperature control range	_	4 to 40°C, 1°C steps			
Cell temperature setting range	_	(Room temperature – 10°C) to 40°C			
Functions	Four-wavelength detection, wavelength scanning				
Safety measures	Liquid-leakage sensor				
Operating temperature range	4 to 35°C				
Dimensions and weight	W260 × D420 × H210 mm, 16 kg	W260 × D420 × H210 mm, 18 kg			

ELSD-LT III



	ELSD-LT III (S228-65900-58)
Nebulizing Method	Siphon splitting
Light Source	Semiconductor laser
Detector	Photodiode
Temperature Setting Range	Room temperature to 100 °C
Nebulizer Gas	Air or nitrogen*
Mobile Phase Flow Rate (Standard Nebulizer)	0.2 to 2 mL/min
Operating Temperature Range	4 to 35 °C
Operating Humidity Range	20 to 85 %
Dimensions and weight	W 250 × D 530 × H 330 mm, 15.5 kg

^{*} Supply gas at a pressure of about 350 kPa. An air compressor may also be used.

A filter (P/N: S228-45528-92) is also available for filtering out moisture and other matter from the compressor.

Main Optional Accessories

Solvent Delivery Units

Part Name	P/N	Description
FCV-11AL	S228-65611-58	This is the mobile phase selection valve (3 flow lines). An FCV-11AL connection kit is required to connect to an FCV-11AL unit.
FCV-11ALS	S228-65610-58	This is the mobile phase selection valve (1 flow line). An FCV-11AL connection kit is required to connect to an FCV-11AL unit.
FCV-11AL Connection Kit	S228-56249-41	This kit includes connector cables and other items necessary for connecting FCV-11AL and FCV-11ALS units.
780 μL Mixer Kit	S228-57313-41	This parts set includes a mixer and tubing for using TFA or other UV-absorbing substance as a mobile phase.
2 mL Mixer Kit	S228-57313-42	This parts set includes a mixer and tubing for using TFA or other UV-absorbing substance as a mobile phase.
Compatible Volume System Kit	S228-57796-42	This kit decreases the system volume to 650 µL.
Low Volume System Kit	S228-57796-43	This kit decreases the system volume to 460 µL.

Autosamplers

Part Name	P/N	Description
50 μL Sample Loop	S228-56074-44	This sample loop is used for injecting 50 µL volumes. (Standard configuration parts of LC-2060)
100 µL Sample Loop	S228-56074-42	This sample loop is used for injecting 100 µL volumes. (Standard configuration parts of LC-2050)
Optional 500 µL Sample Loop	S228-45405-41	This increases the injection volume to 500 µL.
Optional 2 mL Sample Loop	S228-45405-42	This increases the injection volume to 2 mL.
UHPLC Fitting (set of 1)	S228-56867-41	Fitting for inlet to high-pressure capacity column
UHPLC Fitting (set of 10)	S228-56867-43	Fitting for inlet to high-pressure capacity column
Sample Rack	S228-55735-41	Additional sample rack
Plate for 1 mL Sample Vials (set of 2)	S228-56197-41	Plate used to place 84 1 mL sample vials
Plate for 1.5 mL Sample Vials (set of 2)	S228-50830-92	Plate used to place 54 1.5 mL sample vials
Plate for 4 mL Sample Vials (set of 2)	S228-56197-42	Plate used to place 28 4 mL sample vials
Metal plate for 1.5ml Sample Vials (set of 1)	S228-61615-42	Plate used to place 54 1.5 mL sample vials

Column Ovens

Part Name	P/N	Description
Column Clamp ASSY B5	S228-15617-91	This set of clamps is for adding a column with an outside diameter between 6.4 and 9.5 mm.
Column Clamp ASSY B8	S228-15617-92	This set of clamps is for adding a column with an outside diameter between 9.5 and 12.7 mm.
FCV-14AH	S228-65614-58	Automatic column switching valve with 6 positions and 7 ports which is usable at a pressure of 34.3 MPa max.
FCV-34AH	S228-45185-41	Automatic column switching valve with 6 positions and 7 ports which is usable at a pressure of 100 MPa max.
FCV Mounting Kit	S228-55765-42	This parts kit is used to secure an FCV-14AH/ 34AH unit inside the column oven.
CMD	S228-37281-41	This column management device is used to record information about columns.
CMD Cable	S228-39991	This cable is used to connect between the CMD and main units.

UV Detectors

Part Name	P/N	Description
Recycle Valve	S228-56808-41	This low-pressure flow-line selection valve is used to recycle mobile phase.
Flow Cell for UV Detectors	S228-56167-41	This cell is compatible with conventional analysis. (Standard configuration parts of LC-2050)
UHPLC Cell for UV Detectors	S228-45621-41	This cell is compatible with UHPLC analysis. (Standard configuration parts of LC-2060)
Semi-Micro Cell for UV Detectors	S228-45605-46	This cell is compatible with semi-micro analysis.

PDA Detectors

Part Name	P/N	Description
W Lamp ASSY for PDA Detectors	S228-57110-41	This assembly includes a tungsten lamp and its socket used for high-sensitivity analysis in the long-wavelength region.
Flow Cell for PDA Detectors	S228-42593-43	This cell is compatible with conventional analysis. (Standard configuration parts of LC-2050)
High-Speed Cell for PDA Detectors	S228-45618-54	This cell is compatible with fast analysis. (Standard configuration parts of LC-2060)
Semi-Micro Cell for PDA Detectors	S228-45605-47	This cell is compatible with semi-micro analysis.

Other Options

Part Name	P/N	Description
Earthquake Reinforcement Kit	S228-56298-41	This kit is used to reinforce how the reservoir tray is attached.
1 L Mobile Phase Bottles (set of 5)	S228-38583-42	This is a set of five one-liter reservoir bottles for holding mobile phases.
Optional Detector Attachment Kit	S228-56245-41	This kit contains a top plate and reservoir tray for installing an additional detector.
Optional Optical Board	S228-55518-41	This board is used to install additional connectors for optical link cables. It is used to install fluorescence detector RF-20A series and other detectors.
Camera ASSY for Autosampler	S228-55517-41	This camera is installed inside autosamplers. It allows you to monitor the needle action via the computer screen.
Optional AD Board	S228-55519-41	This is an analog-digital converter board. It is used to input the detector signal as an analog signal, such as when a non-Shimadzu detector is connected.
Touch Panel Protecting Sheet	S228-59212-41	Protecting sheet for touch panel.
Upgrade Kit UV	S228-58993-41	Kit for upgrade from LC-2050 (UV model with sample cooler) to LC-2060.
Upgrade Kit PDA	S228-58993-42	Kit for upgrade from LC-2050 (PDA model) to LC-2060.
Smart Automation Kit (4-mobile phase)	S228-26004-44	This kit includes FCV-14AH for up to six columns switching and other parts.
Smart Automation Kit (7-mobile phase)	S228-26004-43	This kit includes FCV-14AH for up to six columns switching, FCV-11AL for solvent delivery of seven mobile phases and other parts.



- Automated support functions utilizing digital technologies, such as M2M, IoT, and Artificial Intelligence (AI), that enable higher productivity and maximum reliability.
- Allows a system to monitor and diagnose itself, handle any issues during data acquisition without user input, and automatically behave as if it were operated by an expert.
- Supports the acquisition of high quality, reproducible data regardless of an operator's skill level for both routine and demanding applications.

Analytical Intelligence logo, LabSolutions and LCsolution are trademarks of Shimadzu Corporation. Hastelloy is a registered trademark of Haynes International, Inc.



Shimadzu Corporation www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.