

Specifications	ID Number	Net Price	QTY	Total Net
<p>Agilent chiller:</p> <ul style="list-style-type: none"> - Electrical Requirements: 220-240VAC, 50Hz, 1450 VA - Voltage $\pm 10\%$, Frequency $\pm 1\text{Hz}$ - Cooling capacity: - Axial system: 1 kW - Radial system: 200W - Flow rate: 1.1 L/min (0.3 gpm) minimum - Recommended inlet temperature: 20°C (68°F) <ul style="list-style-type: none"> Maximum for Agilent 720 Series ICP-OES: 30°C (86°F) Maximum for Agilent 710 Series ICP-OES: 4°C (113°F) - Minimum inlet pressure: 55 kPa (8.0 psi) - Maximum inlet pressure: 310 kPa (45 psi) - Weight: 181 lb. (82 kg) - Dimensions: 368 x 575 x 702 mm Working Temperature Range °C -10° to 40° 	G8481A			

www.aratahiz.com