

# NEW DVE™ Viscometer

our most affordable digital viscometer

## New User Interface

- Keypad control
- Sharp viewing screen for close up or distance viewing

## No Calculations Required

- Direct reading of viscosity

## Displayed Info:

- Viscosity (cP, P, mPa•s or Pa•s)
- % Torque
- Speed/Spindle

## Easy-to-Use

### Range:

Push for determining full scale range (FSR) viscosity

### 18 Speeds

for complete range capability

### Bubble Level

conveniently located for easy adjustment

**Accuracy:  $\pm 1.0\%$  of range**

**Repeatability:  $\pm 0.2\%$**



## What's Included?

Instrument  
6 spindles (RV/HA/HB)  
or 4 spindles (LV)  
Spindle Guard Leg\*  
Lab Stand (Model A)  
Carrying Case

\*Not applicable to HA or HB torque models

## Optional Accessories

Viscosity Standards  
RV/HA/HB-1 Spindle  
Quick Action Lab Stand  
Temperature Bath  
Small Sample Adapter  
UL Adapter  
Thermosel  
DIN Adapter  
Quick Connect/Extension Links

| MODEL | VISCOSITY RANGE<br>cP(mPa•s) |      | SPEEDS |                      |
|-------|------------------------------|------|--------|----------------------|
|       | Min.                         | Max. | RPM    | Number of Increments |
| DVELV | 1†                           | 2M   | .3-100 | 18                   |
| DVERV | 100††                        | 13M  | .3-100 | 18                   |
| DVEHA | 200††                        | 26M  | .3-100 | 18                   |
| DVEHB | 800††                        | 104M | .3-100 | 18                   |

† 1 cP achieved with UL Adapter accessory. 15 cP on LV with standard spindles.

†† Minimum viscosity is achieved with optional RV/HA/HB-1 spindle.

M = 1 million cP = Centipoise mPa•s = Millipascal•seconds