



- LED display of Speed & Time
- Digital control of parameters
- Automatic soft start
- Ten Deceleration rates
- One Memory / Preset Program
- 1-90 second Pulse (Quick spin)
- CE Marked (see safety features)
- Sensors Tacho, Imbalance & Lid
- Stainless steel bowl & barrier ring

• **Complete with Rotor & Accessories**



Speed 1000-4000 Rpm **Max Rcf** 1400G **Timer** 0-30mins **Size** 23x26x39Cm **Weight** 13.75Kg

1030 Rotors



4444

Double sample block and holder



4448

Bi Cyto Cell (alternative)

