

SK-O330-Pro/SK-L330-Pro Orbital and Linear Shakers



*Shown with 1191821 universal attachment sold separately (see overleaf)

The Camlab Choice SK-0330/SK-L330-Pro Orbital and Linear shakers provide smooth and quiet orbital or linear mixing for a wide range of solutions, culture tubes and flasks. They are compact in size and have an electronic timer for consistent shaking between batches. The shaker also has a continuous run mode. Free software ShakerPC is available to control the unit through an RS232 interface.


Specifications

Part Number	Linear Shaker	SK-L330-Pro	1189971
	Orbital Shaker	SK-O330-Pro	1189973
Speed range	Linear Shaker	SK-L330-Pro	100-350 rpm
	Orbital Shaker	SK-O330-Pro	100-500 rpm
Orbital diameter			10mm
Frequency			50/60Hz
Max load			7.5Kg
Dimensions			420 x 370 x 100
Weight			13.5Kg
Protection class			IP21
Power			30W



Features

- Compact size
- LCD display
- Electronic Timer 1-1199 mins
- Free ShakerPC software



camlab

www.camlab.co.uk

Accessories

Cat No

1191821

Universal attachment with 4 bars for SKO330/SKL330 Orbital/Linear Shaker

