



Be Right™



## Manual Titrators

Get accurate ( $\pm 1\%$ ), convenient manual titrations performed without the bulk, fragility, or waste of a conventional glass burette, using the Hach Digital Titrator. Test for eighteen parameters including chlorine, hardness, and iron. Uses interchangeable titrant cartridges.

Part Number	Reagents included	Includes	Carrying case included	Accessories Included	USD Price
1690001	No	Digital Titrator only	No	No	Contact Hach
1690002	No	Digital Titrator device only	Yes	No	Contact Hach
2270900	No	Digital Titrator with accessories in carrying case	Yes	Yes	Contact Hach