



WIZARD

Vortex Mixer

Exclusive and compact vortex mixer with settable speed, provided with the revolutionary IR sensor or continuous mode



IR Sensor, New Conception of Vortex Mixers

Vortex mixer featuring the highly innovative and patented IR technology to reduce repetitive strains. IR sensor and Continuous mode available. Large choice of complementary accessories to accommodate tubes, plates and flasks.

- Vibration starts when tube is detected by IR (VELP patent)
- IR sensor and continuous modes
- Variable speed operation, electronically controlled up to 3000 rpm
- Optimum resistance to chemicals from technopolymer casing
- Excellent stability with zinc alloy base
- IP 42 protection against particles and liquids
- No fatigue as no pressure is required
- Eliminates repetitive strain injury



INSTRUMENT - CODE

WIZARD 100-240 V / 50-60 Hz **F202A0175**
 * Included: Mixing cup head

TECHNICAL DATA

CONSTRUCTION MATERIAL	Zinc alloy and technopolymer
SUPPORT SYSTEM	3 anti-sliding feet
OPERATIONAL MODE	Sensor / Continuous
SPEED SETTING	Analog
SPEED CONTROL	Electronic
STIRRING SPEED	0 - 3000 rpm
ORBITAL DIAMETER	4,5 mm
TYPE OF MOVEMENT	Orbital
ELECTRONIC PROTECTION DEGREE CEI EN 60529	IP 42
POWER	15 W
WEIGHT	2.2 Kg (4.9 lb)
DIMENSIONS (WxHxD)	180x70x220 mm (7.1x2.8x8.7 in)



We reserve the right to make technical alterations
 We do not assume liability for errors in printing, typing or transmission

VELP Official Partner

DESIGNED AND MANUFACTURED IN ITALY



ITALY - HQ
 Via Stazione 16
 20865 Usmate (MB) Italy
 Tel. +39 039 628811
 velpitalia@velp.com

INDIA
 velpindia@velp.com

USA
 155 Keyland Court, Bohemia
 NY 11716 - U.S.
 Tel. +1 631 573 6002
 velpusa@velp.com

CHINA
 Xinlong Rd Building 28, Lane 1333
 Shanghai city - China
 Tel. +86 180 1755 7329
 velpchina@velp.com