

Specifications

Model	ET-24-10 LN2	ET-24-10 Doble
Part Number	21127411	20938931

Maximum Storage Capacity

Number of Canisters	10	10
Number of 1/4 cc straws 24/1.2 ml cane	–	–
Number of 1/4 cc straws 1 Level Bulk	8190	8190
Number of 1/4 cc straws 2 Level Bulk	–	–

Performance

Ln2 Capacity L	26	8
Static Evaporation Rate L/day*	0.35	0.35
Static Holding Time days	74	23
Normal Working Duration days	44	18

Unit Dimensions

Neck Opening mm	97	97
Overall Height mm	480	480
Outer Diameter mm	473	500
Canister Height mm	127	127
Canister Diameter mm	67	67
Weight Empty kg	16.5	16.5
Weight Full kg	43.2	45.7

* Static evaporation rate and static holding time are nominal. Actual rate and holding time will be affected by the nature of the container use, atmospheric conditions, and manufacturing tolerances.

TWO Year Vacuum Warranty from date of invoice