

Agilent 1100 G1330B ALS Therm



[Lease/Rent](#)

[Refer to a Friend](#)

Agilent 1100 G1330B ALS Therm

Listing #: LS-005003

Category:

Chromatography, Liquid Chromatography (HPLC)

Brand: Agilent

Model: G1330B

Availability: In Stock

Condition: Refurbished

Specifications

The Agilent 1100 G1330B ALS Thermostat is designed to provide precise temperature control for autosampler sample trays in Agilent 1100 Series HPLC systems. Integrated directly with Agilent 1100 Series autosamplers, the G1330B supports temperature-sensitive analyses by reducing variability caused by ambient temperature fluctuations. This module is commonly used in pharmaceutical analysis, chemical research, food and beverage testing, and quality control laboratories where reliable sample conditioning is required.

This Agilent 1100 G1330B ALS Thermostat is offered fully refurbished and tested, in working condition, making it suitable for immediate laboratory use.

Includes:

- Agilent 1100 G1330B ALS Thermostat
- Power Cable
- 30-Day Warranty Coverage

Disclaimer: All brand names and trademarks are the property of their respective owners and are used for identification only. We are not an authorized dealer. Software, if included, is not transferable or upgradable.